



REACHING YOUR LAB
IS OUR MISSION



ASIAN SCIENTIFIC INSTRUMENTS

Asian Estate : 4-1-10, Tilak Road Opposite Posnett Bhavan Hyderabad, Telangana-500 001

Ph. : 93930-01939, 90522-23905, 98480-31317, 93472-45096

✉ info@asianhyd.com

CONTENTS

A

| | |
|--------------------------------|----|
| Analytical Funnels | 24 |
| Animal Cage | 19 |
| Animal Cage Twin Grill | 19 |
| Aspirator Bottles | 13 |
| Atomic Model Set | 27 |
| Atomic Model Set (Euro Design) | 27 |

B

| | |
|------------------------------|----|
| Beakers | 1 |
| Beakers (Euro Design) | 1 |
| Beakers (Printed Graduation) | 1 |
| Beakers (Stackable) | 2 |
| Buchner Funnels | 24 |
| Burette | 2 |

C

| | |
|--------------------------------|----|
| Cap for Test Tube | 47 |
| Carboy | 12 |
| Carboy with Stop Cock | 12 |
| Carrier Tray | 69 |
| Centrifuge Tube Box | 45 |
| Centrifuge Tube Conical Bottom | 44 |
| Centrifuge Tube Round Bottom | 46 |
| Collection Vial | 48 |
| Conical Flask | 14 |
| Conical Measures | 64 |
| Coplin Jar | 38 |
| Coplin Jar (Euro Design) | 38 |
| Cross Connectors | 17 |
| Cryo Box (P.C) | 50 |
| Cryo Box (P.P) | 50 |
| Cryo Rack | 49 |
| Cryo Vial | 49 |
| Crystal Model Set | 37 |
| Culture Plate with Cover | 26 |

D

| | |
|---------------------------------|----|
| Desiccator (Plain) | 21 |
| Desiccator (Plain) - All Clear | 21 |
| Desiccator (Vaccum) | 20 |
| Desiccator (Vaccum) - All Clear | 21 |
| Draining Basket | 59 |
| Draining Rack | 70 |
| Draining Tray | 59 |
| Dropping Bottles | 7 |
| Dropping Bottles (Euro Design) | 7 |

F

| | |
|------------------------------|----|
| Fast Releaser Pipette Pump | 6 |
| Fisher Clamp / Burette Clamp | 15 |
| Flask Stand | 16 |
| Float Rack | 41 |
| Forceps | 66 |
| Funnel Holder | 15 |
| Funnel Short Stem | 25 |
| Funnels-Long Stem | 24 |

H

| | |
|--------------------------|----|
| Heavy Duty Vacuum Bottle | 12 |
|--------------------------|----|

I

| | |
|-----------------------------|----|
| Industrial Funnels | 25 |
| Instrument Sterilizing Tray | 69 |

K

| | |
|-----------------|----|
| Kipps Apparatus | 23 |
|-----------------|----|

L

| | |
|--------------------|----|
| L Connectors | 17 |
| L Shaped Spreader | 67 |
| Laboratory Tray | 68 |
| Leclanche Cell Pot | 72 |
| Long Handle Scoop | 65 |

M

| | |
|-------------------------------------|----|
| Magenta Box | 70 |
| MCT Box | 44 |
| MCT Rack | 42 |
| MCT Twin Rack | 42 |
| Measuring Cylinders Hexagonal | 63 |
| Measuring Jugs | 63 |
| Measuring Jugs (Euro Design) | 63 |
| Measuring Jugs (Printed Graduation) | 64 |
| Medicine Cup | 72 |
| Micro Pestle | 67 |

M

| | |
|---|----|
| Micro Test Plate | 26 |
| Micro Tip Box | 51 |
| Micro Tip Box (Euro Design) | 51 |
| Molecular Model Set - | |
| Advance Level Chemistry Set | 31 |
| Atomic Orbital Set - 14 Model Kit | 37 |
| Bio Chemistry Set - Student | 31 |
| Bio Chemistry Set - Teacher | 30 |
| Buckminster Fullerene | |
| "Bucky Ball" - C60 | 35 |
| Cholesterol (C ₂₇ H ₄₆ O) | 35 |
| Diamond | 33 |
| Fat (Glycerol Tristearate) (C ₅₇ H ₁₁₀ O ₆) | 34 |
| Glucose (C ₆ H ₁₂ O ₆) | 34 |
| Graphite | 31 |
| Inorganic / Organic Set - Student | 29 |
| Inorganic / Organic Set - Teacher | 28 |
| Molecular Geometry Set (VSEPR) | 36 |
| Organic Set - Student | 28 |
| Organic Set - Teacher | 27 |
| Organic Stereo Chemistry Set | |
| - Student | 30 |
| Organic Stereo Chemistry Set | |
| - Teacher | 29 |
| Shapes of Molecules - 8 Model Kit | 36 |
| Sodium Chloride - Compact 27 Atoms | 32 |
| Sodium Chloride - Open 27 Atoms | 32 |
| Sodium Chloride - Open 36 Atoms | 32 |
| Scroce (C ₁₂ H ₂₂ O ₁₁) | 34 |
| Sulphur | 35 |
| Water/Ice (H ₂ O) | 33 |
| Zinc Sulphide | 33 |
| Multi Purpose Rectangular Storage | |
| Box | 71 |
| N | |
| Narrow Mouth Bottle | 11 |
| Nestler Cylinder Stand | 54 |
| O | |
| OAK Ridge Centrifuge Tubes | 46 |

CONTENTS

P

| | |
|---|----|
| Pasteur Pipette | 7 |
| PCR Tube Rack | 43 |
| PCR Tube Rack with Hinges- Autoclavable | 43 |
| Perforated Polypropylene Disc for Desiccator | 22 |
| Petri Dish | 26 |
| Petri Dish (Culture) | 26 |
| Petri Dishes Disposable Aseptic | 52 |
| Petri Dishes Disposable ETO Sterile | 52 |
| Petri Dishes Disposable Gamma Sterile | 52 |
| Pinch Clip | 18 |
| Pipette Box | 61 |
| Pipette Pump | 6 |
| Pipette Stand (94 Pipette-Rotary) | 60 |
| Pipette Stand (Horizontal) | 60 |
| Pipette Stand (Vertical) | 60 |
| Pipette Tray | 61 |
| Pipette Tray with 4 Compartments | 61 |
| Plantation Pots | 68 |
| Plastic Stopper | 47 |
| Pneumatic Trough | 68 |
| Policemen Stirring Rod | 67 |
| Powder Funnels | 25 |

R

| | |
|---|----|
| Rack for Centrifuge Tube (Folding)- Space Saver | 45 |
| Rack for Micro Centrifuge Tube (Folding) | 58 |
| Rack for Micro Centrifuge Tube | 55 |
| Rack for Petri Dishes | 52 |
| Rack for Scintillation Vial | 48 |
| Reagent Bottles (Narrow Mouth)-Amber | 8 |
| Reagent Bottles (Narrow Mouth) - HDPE | 8 |
| Reagent Bottles (Narrow Mouth) - LDPE | 9 |
| Reagent Bottles (Narrow Mouth) - PP | 8 |
| Reagent Bottles (Narrow Mouth) - (Printed Graduation) - PP | 9 |
| Reagent Bottles (Wide Mouth)-Amber | 10 |
| Reagent Bottles (Wide Mouth) - HDPE | 10 |
| Reagent Bottles (Wide Mouth) - LDPE | 10 |
| Reagent Bottles (Wide Mouth) - PP | 9 |
| Reagent Bottles (Wide Mouth) - (Printed Graduation) - PP | 11 |
| Reagent Reservoir | 62 |
| Retort Stand | 20 |
| Rotatable Multi Rack | 53 |
| Rotatable Multi Rack-Expandable | 54 |

S

| | |
|--|----|
| Safety Goggles with Vent | 66 |
| Sample Container (Press & Fit Type) | 72 |
| Scintillation Vial | 47 |
| Scoop | 65 |
| Separatory Funnel | 64 |
| Separatory Funnel Holder | 15 |
| Silicon "O" Ring for Desiccator | 22 |
| Simple Cell Pot | 71 |
| Slide Box | 39 |
| Slide Box with Hinges | 40 |
| Slide Draining Tray | 41 |
| Slide File | 40 |
| Slide Mailer | 38 |
| Slide Mailer (Transparent) | 39 |
| Slide Storage Rack | 41 |
| Slide Tray | 40 |
| Soft Loops | 67 |
| Spatula | 65 |
| Specimen Jar (Gas Jar) | 14 |
| Staining Box | 62 |
| Stirrer | 66 |
| Stop Cock For Aspirator Bottle & Carboy | 13 |
| Stop Cock For Desiccator | 22 |
| Stop Cock For Kipps Apparatus | 23 |
| Stop Cocks | 23 |
| Storage Boxes | 71 |
| Storage Vial | 47 |
| Straight Connectors (Tapered) for Tubbing | 17 |

T

| | |
|---------------------------------------|----|
| T Connectors | 16 |
| Test Tube Baskets | 59 |
| Test Tube Baskets (with Cover) | 59 |
| Test Tube Peg Rack | 57 |
| Test Tube Rack | 57 |
| Test Tube Stand | 54 |
| Test Tube Stand (3 Tier) P.C | 56 |
| Test Tube Stand (3 Tier) P.P | 56 |
| Test Tube Stand (Round) | 55 |
| Test Tube Stand (Wire Pattern) | 58 |
| Test Tube Stand (Wire Pattern-Fix) | 58 |
| Test Tube with Screw Cap | 46 |
| Thisle Funnel Set for Kipps Apparatus | 23 |
| Three Tier Micro Tube Rack | 42 |

U

| | |
|--|----|
| Unequal Straight Connectors Tapered for Tubbing | 18 |
| Unequal Straight Connectors Untapered for Tubbing | 48 |
| Universal Combi Link Rack | 53 |
| Universal Multi Rack | 53 |
| Universal Reagent Reservoir | 62 |
| Utility Tray | 69 |

V

| | |
|------------------|----|
| Volumetric Flask | 51 |
|------------------|----|

W

| | |
|-------------------------|----|
| Wash Bottles | 13 |
| Wash Bottles (New Type) | 14 |
| Water Bottle | 19 |
| Wide Mouth Bottle | 11 |

Y

| | |
|--------------|----|
| Y Connectors | 16 |
|--------------|----|



ASIAN SCIENTIFIC INSTRUMENTS
SUPER MARKET FOR LAB ITEMS

BEAKERS

Material : PP Material Conforming US FDA 21CFR

PolyLab Beakers, moulded in Polypropylene, have excellent clarity and very good chemical resistance. These autoclavable Beakers have prominently raised graduations on the outer surface. Pouring a solution becomes very convenient because of its tapered spout.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 11101 | 50ml | 39269099 |
| 11102 | 100ml | 39269099 |
| 11103 | 250ml | 39269099 |
| 11104 | 500ml | 39269099 |
| 11105 | 1000ml | 39269099 |
| 11106 | 2000ml | 39269099 |



BEAKERS (EURO DESIGN)

Material : PP Material Conforming US FDA 21CFR

PolyLab Beakers, moulded in Polypropylene, have excellent clarity and very good chemical resistance. These autoclavable Beakers have prominently raised graduations on the outer surface. Pouring a solution becomes very convenient because of its tapered spout.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 11110 | 25ml | 39269099 |
| 11111 | 50ml | 39269099 |
| 11112 | 100ml | 39269099 |
| 11113 | 250ml | 39269099 |
| 11114 | 500ml | 39269099 |
| 11115 | 1000ml | 39269099 |
| 11116 | 2000ml | 39269099 |
| 11117 | 5000ml | 39269099 |



BEAKERS (PRINTED GRADUATION)

Material : PP Material Conforming US FDA 21CFR

Like other PolyLab beakers, these beakers also, moulded in polypropylene have excellent clarity as well as very good chemical resistance. The prominent feature of these beakers is the easy to read printed graduations. These beakers being in PP are autoclavable.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 11151 | 50ml | 39269099 |
| 11152 | 100ml | 39269099 |
| 11153 | 250ml | 39269099 |
| 11154 | 500ml | 39269099 |
| 11155 | 1000ml | 39269099 |
| 11156 | 2000ml | 39269099 |
| 11157 | 5000ml | 39269099 |
| 11158 | 10000ml | 39269099 |



SINCE 1966



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

BEAKERS (STACKABLE)

Material : PP Material Conforming to US FDA 21CFR

Adding to the existing range of PolyLab Beakers, these beakers are also molded in Polypropylene, have excellent clarity and very good chemical resistance. The unique feature of these beakers is its tapered stacking style to make efficient use of your storage space.



| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 11161 | 50ml | 39269099 |
| 11162 | 100ml | 39269099 |
| 11163 | 250ml | 39269099 |
| 11164 | 600ml | 39269099 |



BURETTE

Material : Polymethyl Methacrylate/PMP

Polylab Burette has a transparent acrylic pipe and a leakproof PMP stopcock filled with a self lubricating PTFE plug. Each Burette is individually calibrated to give maximum accuracy. Acrylic body is unaffected by dilute mineral acids & bases except Hydrofluoric Acid & Ammonium Hydroxide. The pH range of 1-14 is ideal and can be used up to 5% wt/wt of acids or 30% bases.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 11501 | 25ml | 39269099 |
| 11502 | 50ml | 39269099 |
| 11503 | 100ml | 39269099 |

✉ info@asianhyd.com





ASIAN SCIENTIFIC INSTRUMENTS
SUPER MARKET FOR LAB ITEMS

PIPETTE PUMP

Material : Acrylonitrile Butadiene Styrene/PP/Silicon

Polylab presents a fast release pipetting device for precise pipetting & zip quick emptying. The Knurled thumbwheel can be rotated to draw up or dispense liquids while the releaser may be pushed in for zip-quick emptying. The moulded chuck has a silicon rubber collet inside to hold various pipettes of standard sizes. It resists acids & alkalis and can be easily disassembled for through cleaning.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 20001 | 2ml | 39269099 |
| 20002 | 10ml | 39269099 |
| 20003 | 25ml | 39269099 |



FAST RELEASER PIPETTE PUMP

Material : PP/Silicon

These fully autoclavable fast release pipette provides comfort, quick and accurate pipetting. Its unique design reduces the risk of injury and undue pressure on hand and wrist while repetitive pipetting. These lightweight pipette pump has a fast release trigger that dispenses the entire pipette contents while knurled thumbwheel can be rotated to draw up or dispense liquids. The molded chuck has a silicon rubber collet inside to hold various pipettes of standard sizes. It resists acids and alkalis and can be easily disassembled for cleaning.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 11101 | 50ml | 39269099 |
| 11102 | 100ml | 39269099 |
| 11103 | 250ml | 39269099 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS



PASTEUR PIPETTE

Material : Low Density Polyethylene

These disposable pipettes made in low density polyethylene have a bulb and graduations also. These pipette have excellent reproducibility of the no. of drops per milliliters thus making them a best option for distributing aliquots of liquids portions. The long slender tip guarantees precision transfer of the sample.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 20101 | 1ml | 39269099 |
| 20102 | 3ml | 39269099 |



DROPPING BOTTLES

Material : Low Density Polyethylene

Blown in Low Density Polyethylene, Polylab Dropping Bottles dispense small uniform drops of liquids when pressed lightly.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 22101 | 60ml | 39233090 |
| 22102 | 125ml | 39233090 |

SINCE 1966



DROPPING BOTTLES (EURO DESIGN)

Material : Low Density Polyethylene/PP

These bottles made in low density polyethylene with PP caps. Dispense small uniform drops of variety of liquids. The unique feature of these bottles is that they are available with screw cap with lid. Repeatedly produce a measured and accurate drop when squeezed.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 22200 | 8ml | 39233090 |
| 22201 | 15ml | 39233090 |
| 22202 | 30ml | 39233090 |
| 22203 | 60ml | 39233090 |
| 22204 | 125ml | 39233090 |
| 22205 | 250ml | 39233090 |

✉ info@asianhyd.com





ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

REAGENT BOTTLES (NARROW MOUTH)-PP

Material : PP Material Conforming to US FDA 21CFR.

These Narrow Mouth Reagent Bottles are available in Polypropylene. The PP bottles are rigid, translucent and have excellent contact clarity. These bottles can be autoclaved filled as well as empty. The Polypropylene Cap has a built in seal ring to make these bottles leakproof.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 33251 | 4ml | 39233090 |
| 33252 | 8ml | 39233090 |
| 33253 | 15ml | 39233090 |
| 33254 | 30ml | 39233090 |
| 33301 | 60ml | 39233090 |
| 33302 | 125ml | 39233090 |
| 33303 | 250ml | 39233090 |
| 33304 | 500ml | 39233090 |
| 33305 | 1000ml | 39233090 |



REAGENT BOTTLES (NARROW MOUTH)-HDPE

Material : HDPE Natural/PP

These Narrow Mouth Reagent Bottles are available in high density Polypropylene with PP Caps are highly reliable & durable. Rigid & translucent, these bottles are not autoclavable. The seal ring inside the cap make these bottles totally leak proof.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 33401 | 4ml | 39233090 |
| 33402 | 8ml | 39233090 |
| 33403 | 15ml | 39233090 |
| 33404 | 30ml | 39233090 |
| 33405 | 60ml | 39233090 |
| 33406 | 125ml | 39233090 |
| 33407 | 250ml | 39233090 |
| 33408 | 500ml | 39233090 |
| 33409 | 1000ml | 39233090 |



REAGENT BOTTLES (NARROW MOUTH)-AMBER

Material : HDPE Amber/PP Amber

These Amber Narrow Mouth Reagent Bottles available in HDPE with Amber PP caps are highly recommended for reducing UV light transmission to help protect light sensitive reagents. Rigid & opaque, these bottles are not auto clavable. The seal ring inside the cap make these bottles totally leak proof.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 33421 | 4ml | 39233090 |
| 33422 | 8ml | 39233090 |
| 33423 | 15ml | 39233090 |
| 33424 | 30ml | 39233090 |
| 33425 | 60ml | 39233090 |
| 33426 | 125ml | 39233090 |
| 33427 | 250ml | 39233090 |
| 33428 | 500ml | 39233090 |
| 33429 | 1000ml | 39233090 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

Asian Estate : 4-I-10, Tilak Road Opposite Posnett Bhavan Hyderabad, Telangana-500 001

Ph. : 93930-01939, 90522-23905, 98480-31317, 93472-45096

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS



REAGENT BOTTLES (NARROW MOUTH)-LDPE

Material : LDPE/PP

These Narrow Mouth Bottles available in LDPE with PP caps, have excellent flexibility, chemical resistance & impact resistance. Contents are visible. These bottles are not autoclavable. The seal ring inside the cap make these bottles totally leak proof.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 33431 | 8ml | 39233090 |
| 33432 | 15ml | 39233090 |
| 33433 | 30ml | 39233090 |
| 33434 | 60ml | 39233090 |
| 33435 | 125ml | 39233090 |
| 33436 | 250ml | 39233090 |
| 33437 | 500ml | 39233090 |
| 33438 | 1000ml | 39233090 |



REAGENT BOTTLES (NARROW MOUTH) (PRINTED GRADUATION) - PP

Material : PP Material Conforming to US FDA 21CFR

Adding to the existing features of polypropylene narrow mouth reagent bottle, these bottles comes with printed easy to read approximate volume graduations.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 33481 | 8ml | 39233090 |
| 33482 | 15ml | 39233090 |
| 33483 | 30ml | 39233090 |
| 33484 | 60ml | 39233090 |
| 33485 | 125ml | 39233090 |
| 33486 | 250ml | 39233090 |
| 33487 | 500ml | 39233090 |
| 33488 | 1000ml | 39233090 |

SINCE 1966



REAGENT BOTTLES (WIDE MOUTH)-PP

Material : PP Material Conforming to US FDA 21CFR.

These Wide Mouth Reagent Bottles are available in Polypropylene. The PP bottles are rigid, translucent and have excellent contact clarity. These bottles can be autoclaved filled as well as empty. The Polypropylene Cap has a built in seal ring to make these bottles leakproof.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 33311 | 30ml | 39233090 |
| 33306 | 60ml | 39233090 |
| 33307 | 125ml | 39233090 |
| 33308 | 250ml | 39233090 |
| 33309 | 500ml | 39233090 |
| 33310 | 1000ml | 39233090 |



ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

REAGENT BOTTLES (WIDE MOUTH)-HDPE

Material : HDPE Natural/PP

These Wide Mouth Reagent Bottles are available in high density Polypropylene with PP Caps are highly reliable & durable. Rigid & translucent, these bottles are not autoclavable. The seal ring inside the cap make these bottles totally leak proof.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 33441 | 30ml | 39233090 |
| 33442 | 60ml | 39233090 |
| 33443 | 125ml | 39233090 |
| 33444 | 250ml | 39233090 |
| 33445 | 500ml | 39233090 |
| 33446 | 1000ml | 39233090 |



REAGENT BOTTLES (WIDE MOUTH)-AMBER

Material : HDPE Amber/PP Amber

These wide Mouth Reagent Bottles available in HDPE with PP Amber caps are highly recommended for reducing UV light transmission to help protect light sensitive reagents. Rigid & opaque, these bottles are not autoclavable. The seal ring inside the cap make these bottles totally leak proof.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 33461 | 30ml | 39233090 |
| 33462 | 60ml | 39233090 |
| 33463 | 125ml | 39233090 |
| 33464 | 250ml | 39233090 |
| 33465 | 500ml | 39233090 |
| 33466 | 1000ml | 39233090 |



SINCE 1966

REAGENT BOTTLES (WIDE MOUTH)-LDPE

Material : LDPE/PP

These Wide Mouth bottles available in LDPE with PP caps, have excellent flexibility, chemical resistance & impact resistance. Contents are easily visible. These bottles are not autoclavable. The seal ring inside the cap makes these bottles totally leak proof.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 33471 | 30ml | 39233090 |
| 33472 | 60ml | 39233090 |
| 33473 | 125ml | 39233090 |
| 33474 | 250ml | 39233090 |
| 33475 | 500ml | 39233090 |
| 33476 | 1000ml | 39233090 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS



**REAGENT BOTTLES (WIDE MOUTH)
(PRINTED GRADUATION) - PP**

Material : PP Material Conforming to US FDA 21CFR

Adding to the existing features of polypropylene wide mouth reagent bottle, these bottles come with printed easy to read approximate volume graduations.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 33491 | 30ml | 39233090 |
| 33492 | 60ml | 39233090 |
| 33493 | 125ml | 39233090 |
| 33494 | 250ml | 39233090 |
| 33495 | 500ml | 39233090 |
| 33496 | 1000ml | 39233090 |



NARROW MOUTH BOTTLE)

Material : PP Material Conforming to US FDA 21CFR

This autoclavable bottle, made of polypropylene, can be used for storing & transferring chemicals and other reagents safe & sound.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 33500 | 2000ml | 39233090 |
| 33501 | 4000ml | 39233090 |
| 33502 | 8000ml | 39233090 |

SINCE 1966



WIDE MOUTH BOTTLE)

Material : PP Material Conforming to US FDA 21CFR

This autoclavable bottle, made of polypropylene, can be used for storing & shipping chemicals & reagents. Being wide mouthed, it is always convenient to pour in & take out the reagents. The built in depression on two sides in square bottle provide excellent grip while holding this bottle.

| Product Code | Particulars | HSN Code |
|--------------|---------------|----------|
| 33510 | Round 2000ml | 39233090 |
| 33511 | Square 4000ml | 39233090 |
| 33512 | Round 4000ml | 39233090 |



ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

HEAVY DUTY VACCUM BOTTLE

Material : PP Material Conforming to US FDA 21CFR.

This autoclavable bottle, made of polypropylene, has thick walls, is rugged & can be used for repeated applications of vacuum.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 33521 | 2000ml | 39233090 |



CARBOY

Material : PP Material Conforming to US FDA 21CFR.

Polylab Graduated Carboys, made of polypropylene, are much lighter than glass. These autoclavable, dependable & durable Carboys are provided with wide handles making it convenient while carrying or pouring purposes. Ideal for storing & transporting distilled water, culture media & other laboratory solutions.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 34001 | 10Ltrs. | 39233090 |
| 34002 | 20Ltrs. | 39233090 |



SINCE 1966

CARBOY WITH STOP COCK

Material : PP Material Conforming to US FDA 21CFR.

These polypropylene graduated Carboys are autoclavable, dependable & durable. These are also equipped with a leak proof spigot screwed to the body. These carboys are an excellent replacement for glass carboys and are used for storing & dispensing distilled water & other laboratory solutions. Threads should be totally disengaged before autoclaving.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 34051 | 10Ltrs. | 39233090 |
| 34052 | 20Ltrs. | 39233090 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS



ASPIRATOR BOTTLES

Material : PP Material Conforming to US FDA 21CFR.

Polylab Aspirator Bottles are made of Polypropylene and are therefore much lighter than glass. These are autoclavable, durable & dependable. The leakproof spigot screwed to the bottle body requires only a half turn to open or close. These bottles are used for storing or distilled water & other laboratory solutions. Threads should be totally disengaged before autoclaving.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 34100 | 2Ltrs. | 39233090 |
| 34101 | 5Ltrs. | 39233090 |
| 34102 | 10Ltrs. | 39233090 |
| 34103 | 20Ltrs. | 39233090 |



STOP COCK FOR ASPIRATOR BOTTLE & CARBOY

Material : PP / Silicon / PTFE

These stop cock are suitable for all sizes of carboys with stop cock & aspirator bottle.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| SP1117 | Stop Cock | 39269099 |



WASH BOTTLES

Material : Polyethylene

Polylab Wash Bottles are made of Low Density Polyethylene which gives these bottles a translucent & unbreakable character. Solutions can be poured easily by squeezing the bottle lightly. The cap is fitted with a flexible polyethylene delivery tube which can be aimed wherever required.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 36601 | 125ml | 39233090 |
| 36602 | 250ml | 39233090 |
| 36603 | 500ml | 39233090 |
| 36604 | 1000ml | 39233090 |



ASIAN SCIENTIFIC INSTRUMENTS SUPER MARKET FOR LAB ITEMS

WASH BOTTLES (NEW TYPE)

Material : Low Density Polyethylene/PP

These Bottles are also made of Low Density Polyethylene provided with PP caps are translucent & unbreakable. The distinguishing feature that these bottles contains is the built in delivery tube making it a one-piece Wash Bottle. These Bottles are graduated on the outer surface to give an approximate idea about the volume of solution filled or dispensed.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 36607 | 125ml | 39233090 |
| 36605 | 250ml | 39233090 |
| 36606 | 500ml | 39233090 |
| 36608 | 750ml | 39233090 |
| 36609 | 1000ml | 39233090 |



CONICAL FLASK

Material : Polypropylene

Made of Polypropylene, these cone shaped flasks are rigid, translucent & autoclavable. The screw cap keeps the contents of the flask contamination free and make the flask leakproof. Threads should be disengaged before autoclaving.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 38101 | 100ml | 39233090 |
| 38102 | 250ml | 39233090 |
| 38103 | 500ml | 39233090 |



SINCE 1966

SPECIMEN JAR (GAS JAR)

Material : Polystyrene

These transparent Specimen Jars, used for preserving different specimens in the lab are provided with a cap that fits in tight. These jars provided excellent protection to species preserved in it.

| Product Code | Particulars | HSN Code |
|--------------|-----------------|----------|
| 40001 | o 50mm x h150mm | 39231090 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS



FISHER CLAMP / BURETTE CLAMP

Material : Polypropylene

These Clamps are moulded in Polypropylene making it completely corrosion-proof. The spring loaded jack in the middle provides a vice like grip to the Burette and prevents it from slipping. Two tapered wedges are provided to ensure firm grip on rods having dia 1/2" or 3/8".

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 45101 | Singal | 39269099 |
| 45102 | Double | 39269099 |



FUNNEL HOLDER

Material : Polypropylene

Made of Polypropylene, these corrosion free Funnel Holders can hold funnels with dia 3" to 6". The extra plate provided can be placed to hold funnels bearing dia 1" to 3". Two tapered wedges are provided to ensure firm grip on Roads having dia 1/2" or 3/8".

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 45103 | Singal | 39269099 |
| 45106 | Double | 39269099 |

SINCE 1966



SEPARATORY FUNNEL HOLDER

Material : Polypropylene

This non-corrosive Separatory Funnel Holder, moulded in Polypropylene, has a front opening which provides clear view of the solution and allows unobstructed placement of separating funnels. Two tapered wedges are provided to ensure firm grip on rods having dia 1/2" or 3/8".

| Product Code | Particulars | HSN Code |
|--------------|--------------------------|----------|
| 45104 | Separatory Funnel Holder | 39269099 |



ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

FLASK STAND

Material : Polypropylene

Moulded in Polypropylene, this Flask Stand gives a stair-like top view. Every step provides excellent support to the round bottom flask at every point of Contact. This Flask Stand can hold round bottom Flask up to 10 Ltr. capacity. This Stack able Flask Stand can also be steam autoclaved.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 45105 | Flask Stand | 39269099 |



T CONNECTORS

Material : PP Material Conforming to US FDA 21CFR.

These connectors add to the existing range of connectors. The crest & valley structure provided at all the ends of the connectors give excellent leakproof grip to the tube.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 46111 | T 6mm | 39269099 |
| 46113 | T 8mm | 39269099 |
| 46115 | T 10mm | 39269099 |



SINCE 1966

Y CONNECTORS

Material : PP Material Conforming to US FDA 21CFR.

These connectors add to the existing range of connectors. The crest & valley structure provided at all the ends of the connectors give excellent leakproof grip to the tube.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 46112 | Y 6mm | 39269099 |
| 46114 | Y 8mm | 39269099 |
| 46116 | Y 10mm | 39269099 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS



CROSS CONNECTORS

Material : PP Material Conforming to US FDA 21CFR.

These four way connectors prove to be an excellent option for achieving leakproof joints in tubing connections.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 46121 | 6mm | 39269099 |
| 46122 | 8mm | 39269099 |
| 46123 | 10mm | 39269099 |



L CONNECTORS

Material : PP Material Conforming to US FDA 21CFR.

These L shaped connectors also add to the range of different connectors committed to provide excellent leakproof tube arrangements.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 46131 | L 6mm | 39269099 |
| 46132 | L 8mm | 39269099 |
| 46133 | L 10mm | 39269099 |

SINCE 1966



STRAIGHT CONNECTORS (TAPERED) FOR TUBBING

Material : PP Material Conforming to US FDA 21CFR.

Straight connectors with tapered conical shape, with same fitting diameters on both the sides allow connecting tubing having same inner diameters.

| Product Code | Particulars | HSN Code |
|--------------|-----------------------|----------|
| 46141-S1 | with an ID of 3-5mm | 39269099 |
| 46141-S2 | with an ID of 5-7mm | 39269099 |
| 46141-S3 | with an ID of 7-10mm | 39269099 |
| 46141-S4 | with an ID of 9-12mm | 39269099 |
| 46141-S5 | with an ID of 11-14mm | 39269099 |
| 46141-S6 | with an ID of 13-16mm | 39269099 |



ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

UNEQUAL STRAIGHT CONNECTORS TAPERED FOR TUBING

Material : PP Material Conforming to US FDA 21CFR.

Unequal straight connectors with tapered conical shape, with different fitting diameters on both the sides allow connecting tubing having same inner diameters.

| Product Code | Particulars | HSN Code |
|--------------|----------------------------|----------|
| 46151-R1 | with an ID of 4-8/8-12mm | 39269099 |
| 46151-R2 | with an ID of 4-8/12-16mm | 39269099 |
| 46151-R3 | with an ID of 8-12/12-16mm | 39269099 |



UNEQUAL STRAIGHT CONNECTORS UTAPERED FOR TUBING

Material : PP Material Conforming to US FDA 21CFR.

Unequal straight connectors with untapered cylindrical shape, with different fitting diameters on both the sides allow connecting tubing having different inner diameters.

| Product Code | Particulars | HSN Code |
|--------------|----------------------|----------|
| 46151-R1 | with an ID of 7-13mm | 39269099 |



SINCE 1966

PINCH CLIP

Material : PP Material Conforming to US FDA 21CFR.

PolyLab Pinch clip moulded in polypropylene is the best in class ON/OFF ratchet clamp that can be used with tubing having O.D.5mm to 9.5mm. These Pinch clips can be autoclaved.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 46201 | Pinch Clip | 39269099 |



✉ info@asianhyd.com


ASIAN SCIENTIFIC INSTRUMENTS

Asian Estate : 4-I-10, Tilak Road Opposite Posnett Bhavan Hyderabad, Telangana-500 001

Ph. : 93930-01939, 90522-23905, 98480-31317, 93472-45096

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS




ANIMAL CAGE

Material : Polypropylene / Stainless Steel

Available in two different sizes, these caps provide excellent shelter for different animal species in the laboratory. The metallic top provides proper ventilation and has provisions for feeding the species.

| Product Code | Particulars | HSN Code |
|--------------|----------------|----------|
| 47001 | 290x220x140mm. | 39269099 |
| 47002 | 290x220x140mm. | 39269099 |



ANIMAL CAGE TWIN GRILL

Material : Polypropylene / Stainless Steel

Adding in line with the existing range, these Animal Cages are provided with a Stainless Steel grill at the bottom of the cage to provide a clean surface for the species to move on.

| Product Code | Particulars | HSN Code |
|--------------|----------------|----------|
| 48001 | 290x220x140mm. | 39269099 |
| 48002 | 290x220x140mm. | 39269099 |

SINCE 1966



WATER BOTTLE

Material : Polypropylene / Stainless Steel

These water bottles prove to be an excellent accessory for the species kept in the cages.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 47003 | 125 ml | 39233090 |
| 47004 | 250 ml | 39233090 |



ASIAN SCIENTIFIC INSTRUMENTS SUPER MARKET FOR LAB ITEMS

RETORT STAND

Material : Polystyrene/Polypropylene

These non-corrosive, moulded Retort Stands provide excellent replacement for the metal ones. Even the metallic rod has a polypropylene covering making it non-corrosive. The rod measuring 12.5mm x 750mm can be screwed into the heavy base giving way to a perfectly stable Retort Stand.

| Product Code | Particulars | HSN Code |
|--------------|-----------------------|----------|
| 54101 | 225x150mm Side Hole | 39269099 |
| 54102 | 300x200mm Centre Hole | 39269099 |
| 54103 | 300x200mm Side Hole | 39269099 |



DESICCATOR (VACCUUM)

Material : Polystyrene/Polypropylene

Polylab Desiccators are made of Polypropylene and Polycarbonate. These autoclavable vacuum desiccators can hold vacuum upto 740mm of Hg for 24 hrs. without any greasing. The top dome, moulded in rigid and transparent Polycarbonate, gives a crystal clear view of the desiccant placed inside. The Knurled knob on the top provides easy handling of the dome. The lower part is made of Polypropylene. An internal groove is provided on the flange to hold a silicon rubber O-Ring. The Polypropylene stopcock is fitted with a self lubricating PTFE plug. This plug works three way i.e vacuum creation, shutting off & vacuum releasing.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 55204 | 150 mm | 39269099 |
| 55205 | 200 mm | 39269099 |
| 55206 | 250 mm | 39269099 |
| 55207 | 300 mm | 39269099 |
| 55210 | 380 mm | 39269099 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS



DESICCATOR (VACCUM) ALL CLEAR

Material : Polypropylene/Polycarbonate

These Vaccum Desiccators are provided with a transparent polycarbonate bottom that makes it more robust and provides excellent view of the desiccator. All the other features are same as above vaccum desiccators.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 55214 | 150 mm | 39269099 |
| 55215 | 200 mm | 39269099 |
| 55216 | 250 mm | 39269099 |
| 55217 | 300 mm | 39269099 |
| 55220 | 380 mm | 39269099 |



DESICCATOR (PLAIN)

Material : Polypropylene/Polycarbonate

Plain Desiccators are physically the same as Vaccum Desiccator. Only difference is the absence of stop cock. Grease should be applied on the flanges to make it air tight.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 55201 | 150 mm | 39269099 |
| 55202 | 200 mm | 39269099 |
| 55203 | 250 mm | 39269099 |
| 55208 | 300 mm | 39269099 |
| 55209 | 380 mm | 39269099 |

SINCE 1966

DESICCATOR (PLAIN) - ALL CLEAR

Material : Polycarbonate

These plain desiccators are provided with a transparent polycarbonate bottom that makes it more robust and provides excellent view of the desiccators. Rest all other features are same as PP/PC plain desiccators.



| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 55211 | 150 mm | 39269099 |
| 55212 | 200 mm | 39269099 |
| 55213 | 250 mm | 39269099 |
| 55218 | 300 mm | 39269099 |
| 55219 | 380 mm | 39269099 |

Both Vaccum & Plain Desiccators are provided with a thick perforated Polypropylene disc with a big hole in the centre for easy lifting. The tiny holes in the disc provide air troughs & support to the dishes. However porcelain discs are advisable in case of incandescent crucibles.



ASIAN SCIENTIFIC INSTRUMENTS
SUPER MARKET FOR LAB ITEMS

STOP COCK FOR DESICCATOR

Material : Polypropylene / PTFE

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 46151-R1 | with an ID | 39269099 |



SILICON "O" RING FOR DESICCATOR

Material : Silicon

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| SP1111 | 150 mm | 39269099 |
| SP1112 | 200 mm | 39269099 |
| SP1113 | 250 mm | 39269099 |
| SP1114 | 300 mm | 39269099 |
| Sp1122 | 380 mm | 39269099 |




SINCE 1966

PERFORATED POLYPROPYLENE DISC FOR DESICCATOR

Material : Polypropylene

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| SP1118 | 150 mm | 39269099 |
| SP1119 | 200 mm | 39269099 |
| SP1120 | 250 mm | 39269099 |
| SP1121 | 300 mm | 39269099 |
| Sp1123 | 380 mm | 39269099 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS



STOP COCKS

Material : Polypropylene / PTFE

These polypropylene Stop Cocks are fitted with a PTFE plug & have serrated tabulations that accept from 6mm to 10mm tubing. These Stop Cocks are absolutely leakproof and are available in two different options i.e. 2 way & 3 way.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 55301 | 2 Way | 39269099 |
| 55302 | 3 Way | 39269099 |



KIPPS APPARATUS

Material : Polypropylene

Polylab Kipp's Apparatus made of Polypropylene has three main portions. The top portion is fitted with a thistle funnel and a bulb structure through which dilute acid can be poured which reacts with a metal sulphide in the base of the Apparatus. H₂S gas then evolved collects in the middle dome from where it can be released through a leakproof stopcock fitted with a self lubricating PTFE plug. The base is provided with a residual outlet fitted with a cap for easy cleaning.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 55501 | 250 ml | 39269099 |
| 55502 | 500 ml | 39269099 |
| 55503 | 1000 ml | 39269099 |

STOP COCK FOR KIPPS APPARATUS

Material : Polypropylene / PTFE

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| SP1115 | Stop Cock | 39269099 |

THISLE FUNNEL SET FOR KIPPS APPARATUS

Material : Polypropylene / LDPE

| Product Code | Particulars | HSN Code |
|--------------|-----------------|----------|
| Sp1116 | This Funnel Set | 39269099 |



ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

FUNNELS - LONG STEM

Material : PP Material Conforming to US FDA 21CFR.

These long stem funnels have smooth walls, 60° cone angle and external ribs to prevent air lock. Molded in Polypropylene, these clear funnels can also be steam autoclaved.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 57104 | 100ml | 39269099 |
| 57107 | 150ml | 39269099 |



BUCHNER FUNNELS

Material : Polypropylene

Moulded in Polypropylene, these two piece funnels are light weight, easy to clean & autoclavable. The top portion can be snap fitted with the lower one forming a vacuum seal which does not require greasing. The base of the top portion has a perforated filter plate.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 57104 | 70mm | 39269099 |
| 57107 | 110mm | 39269099 |
| 57108 | 90mm | 39269099 |



SINCE 1966

ANALYTICAL FUNNELS

Material : PP Material Conforming to US FDA 21CFR

Designed specifically for analytical chemistry, the body of the funnels is at 60° angle. Moulded in PP, these funnels are provided with internal & external ribs to prevent air blockage.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 57199 | 35 mm | 39269099 |
| 57200 | 50 mm | 39269099 |
| 57201 | 62 mm | 39269099 |
| 57202 | 75 mm | 39269099 |
| 57203 | 100 mm | 39269099 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS



FUNNELS SHORT STEM

Material : PP Material Conforming to US FDA 21CFR

These autoclavable short steam funnel have excellent strength and chemical resistance. Designed for use with standard filter paper. The funnel body is at a 60° angle. Internal and external ribs prevent air locks.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 57211 | 55 mm | 39269099 |



POWDER FUNNELS

Material : PP Material Conforming to US FDA 21 CFR.

These polypropylene funnels are best suited for transferring powdered reagents. The unique feature of these funnels is its parallel stem that minimizes bridging of powder. These funnels are steam autoclavable at 121°C

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 57221 | 60 mm | 39269099 |
| 57222 | 80 mm | 39269099 |
| 57223 | 100 mm | 39269099 |
| 57224 | 150 mm | 39269099 |

SINCE 1966



INDUSTRIAL FUNNELS

Material : PP Material Conforming to US FDA 21 CFR.

Industrial Funnel are designed in such a way that they make filling of carboys, drums and barrels very convenient. The handles of these funnels make their handling very easy.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 57231 | 200 mm | 39269099 |
| 57232 | 250 mm | 39269099 |



ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

PETRI DISH

Material : PP Material Conforming to US FDA 21CFR

These Polypropylene Petri Dishes are clear, unbreakable & autoclavable.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 57301 | 50 mm | 39233090 |
| 57302 | 75 mm | 39233090 |
| 57303 | 100 mm | 39233090 |
| 57304 | 125 mm | 39233090 |
| 57305 | 150 mm | 39233090 |



PETRI DISH (CULTURE)

Material : Polystyrene

These PS Petri Dish with fitting proves excellent for culture exp. Made in polystyrene, this polystyrene, this transparent Petri Dish provides excellent inside view.

| Product Code | Particulars | HSN Code |
|--------------|----------------------|----------|
| 57501 | Petri Dish (Culture) | 39233090 |




MICRO TEST PLATE

Material : Polystyrene

Polylab Micro Test Plates are made of optically clear, non toxic polystyrene. These plates with 96 Flat Bottom and U Bottom wells bear alpha numeric index. Uniform placement of all the wells facilities use of multi channel pipettes.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 57551 | Flat Bottom | 39269099 |
| 57552 | U Bottton | 39269099 |



CULTURE PLATE WITH COVER

Material : Polystyrene

These Clear polystyrene culture plate features deep uniform wells and include lids. Each well hasan alpha-nueric marking. These plates are individually wrapped in a plastic bag.

The 24-well plate measures 5" long x 3-1/2" wide x 7/8" high. Flat bottomed wells are approximately 17.5 mm deep & 16mm in diameter.

The 12-well plate measures 5"long x 3-1/2" wide x 7/8 high. Flat-bottomed wells are

| Product Code | Particulars | HSN Code |
|--------------|-----------------------|----------|
| 57601 | 12 well Plate | 39269099 |
| 57602 | 24 well Plate | 39269099 |
| 57611 | 12 well Plate-Sterile | 39269099 |
| 57612 | 24 well Plate-Sterile | 39269099 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS



ATOMIC MODEL SET

Material : PS/PE

This set consists of moulded balls of different colours & sizes along with connecting lugs of different sizes. These sets are very helpful in modelling different organic & inorganic compounds.

SENIOR SET - This set consists of 185 balls & 370 balls of prescribed colour & sizes packed in a moulded box having separate space for every particular colour & size ball.

JUNIOR SET - This set consists of 37 connecting lugs & 75 balls of different colours & sizes packed in a moulded box.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 58001 | Senior Set | 90230010 |
| 58002 | Junior Set | 90230010 |



ATOMIC MODEL SET (EURO DESIGN)

Material : PS/PE

This set consists of moulded balls of different colours & sizes along with connectors of different sizes. The sizes & shapes of the balls & connecting lugs have been designed in accordance with internationally acceptable standards.

STUDENT SET - This set consists of 47 balls of prescribed colour, size & no. of holes and 35 connectors.

TEACHER SET - This set consists of 100 balls of prescribed colour, size & no. of holes and 86 connectors.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 58011 | Student Set | 90230010 |
| 58012 | Teacher Set | 90230010 |

SINCE 1966

**MOLECULAR MODEL SET - ORGANIC SET
TEACHER**

Material : PS/PE

ELEMENTS

- CARBON (II,III,IV,V) - 38
- OXYGEN (II) - 12
- HYDROGEN (I) - 42
- NITROGEN (IV) - 4
- SULPHUR (II,IV) - 4
- HALOGEN (I) - 8
- METAL (II,VI) - 4

LINKS

- ORBITAL 2D-18
- CONNECTORS - 160
- LINK REMOVERTOOL - 1



| Product Code | Particulars | HSN Code |
|--------------|----------------------|----------|
| 58021 | Organic Set- Teacher | 90230010 |

✉ info@asianhyd.com





BEAKERS

PolyLab
PLASTIC LABWARE

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

MOLECULAR MODEL SET - ORGANIC SET STUDENT

Material : PS/PE


ELEMENTS

CARBON (IV) - 14
OXYGEN (II) - 6
HYDROGEN (I) - 22
NITROGEN (I,III,IV) - 6
SULPHUR (IV, VI) - 2
HALOGEN (I) - 4
METAL (II) - 1
PHOSPHORUS (IV) - 1

LINKS

CONNECTORS - 70
LINK REMOVERTOOL - 1

| Product Code | Particulars | HSN Code |
|--------------|----------------------|----------|
| 58021 | Organic Set- Student | 90230010 |



MOLECULAR MODEL SET - INORGANIC/ ORGANIC SET - TEACHER

Material : PS/PE


ELEMENTS

CARBON (IV, V) - 20
OXYGEN (I, II, IV) - 27
HYDROGEN (I, II) - 14
NITROGEN (III, IV) - 10
SULPHUR (II, IV, VI) - 13
PHOSPHORUS (III, IV, V) - 7
HALOGEN (I) - 8
METAL (I, II, III, IV, VI) - 14

LINKS

ORBITAL 2D - 18
CONNECTORS - 95

| Product Code | Particulars | HSN Code |
|--------------|--------------------------------|----------|
| 58023 | Inorganic/Organic Set- Teacher | 90230010 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS


Asian Estate : 4-I-10, Tilak Road Opposite Posnett Bhavan Hyderabad, Telangana-500 001

Ph. : 93930-01939, 90522-23905, 98480-31317, 93472-45096

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

**MOLECULAR MODEL SET - INORGANIC/
ORGANIC SET - STUDENT**



Material : PS/PE

ELEMENTS

CARBON (IV) - 6
HYDROGEN (I) - 15
NITROGEN (III, IV) - 2
BORON (III) - 1
OXYGEN (I, II, IV) - 27
SULPHUR (II, VI) - 2


PHOSPHORUS (III, V) - 2
HALOGEN (I) - 12
METAL (I, II, III, IV, VI) - 8
sp³ (IV) - 1
dsp³ (V) - 1
d²sp³ - 1

LINKS

ORBITAL 2D - 6
CONNECTORS - 42

| Product Code | Particulars | HSN Code |
|--------------|--------------------------------|----------|
| 58024 | Inorganic/Organic Set- Student | 90230010 |

**MOLECULAR MODEL SET - ORGANIC
STEREO CHEMISTRY SET - TEACHER**



Material : PS/PE

ELEMENTS

HYDROGEN (I, II) - 66
CARBON (IV, V) - 44
OXYGEN (I, II, IV) - 20
NITROGEN (IV, V) - 12
PHOSPHORUS (IV, V) - 5

SULPHUR (II, IV) - 12
IODINE (I) - 2
BROMINE (I) - 2
CHLORINE (I) - 8
METAL (I, II, III, IV, VI) - 8

LINKS

ORBITAL 2D - 12
CONNECTORS - 192
LINK REMOVERTOOL - 1

| Product Code | Particulars | HSN Code |
|--------------|--|----------|
| 58025 | Organic Stereo Chemistry Set - Teacher | 90230010 |



ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

MOLECULAR MODEL SET - ORGANIC STEREO CHEMISTRY SET - STUDENT

Material : PS/PE


ELEMENTS

HYDROGEN (I) - 30
CARBON (IV, V) - 22
OXYGEN (III) - 6
NITROGEN (IV) - 4

BROMINE (I) - 2
CHLORINE (I) - 8
IODINE (I) - 2
METAL (V, VI) - 3

LINKS

ORBITAL 2D - 12
CONNECTORS - 112
LINK REMOVER TOOL - 1



| Product Code | Particulars | HSN Code |
|--------------|--|----------|
| 58026 | Organic Stereo Chemistry Set - Student | 90230010 |

MOLECULAR MODEL SET - BIO CHEMISTRY SET - TEACHER

Material : PS/PE

ELEMENTS

CARBON (II, III, IV) - 71
NITROGEN (II, III, IV) - 36
OXYGEN (I, II) - 47
HYDROGEN (I, II) - 120
SULPHUR (II) - 2
PHOSPHORUS (IV) - 6
METAL (IV, VI) - 2

LINKS

CONNECTORS - 270
LINK REMOVER TOOL - 1



| Product Code | Particulars | HSN Code |
|--------------|-----------------------------|----------|
| 58027 | Bio Chemistry Set - Teacher | 90230010 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS



MOLECULAR MODEL SET - BIO CHEMISTRY SET - STUDENT

Material : PS/PE


ELEMENTS

- CARBON (III, IV) - 24
- OXYGEN (I, II) - 13
- HYDROGEN (I, II) - 28
- NITROGEN (II, III, IV) - 11
- SULPHUR (II) - 1
- PHOSPHORUS (IV) - 1

LINKS

CONNECTORS - 60

| Product Code | Particulars | HSN Code |
|--------------|-----------------------------|----------|
| 58028 | Bio Chemistry Set - Student | 90230010 |



MOLECULAR MODEL SET - ADVANCE LEVEL CHEMISTRY SET - TEACHER

Material : PS/PE

ELEMENTS

- CARBON (IV) - 8
- HYDROGEN (I) - 18
- NITROGEN (III, IV) - 4
- BORON (III, IV) - 2
- OXYGEN (II, IV) - 2
- FLUORINE (I) - 6
- CHLORINE (I) - 6
- SULPHUR (II, IV, VI) - 3
- PHOSPHORUS (III, IV, VI) - 4
- METAL (II, III, IV, VI) - 6

LINKS

CONNECTORS - 42
LOBE - 4

| Product Code | Particulars | HSN Code |
|--------------|-----------------------------|----------|
| 58029 | Advance Level Chemistry Set | 90230010 |



MOLECULAR MODEL SET - GRAPHITE

Material : PS/PE

ELEMENTS

- CARBON (V) - 45

LINKS

CONNECTORS - 72

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 58030 | Graphite | 90230010 |



ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

MOLECULAR MODEL SET - SODIUM CHLORIDE - OPEN 27 ATOMS

Material : PS/PE


ELEMENTS

SODIUM (VI) - 13
CHLORINE (VI) - 14

LINKS

CONNECTORS - 60

| Product Code | Particulars | HSN Code |
|--------------|---------------------------------|----------|
| 58031 | Sodium Chloride - Open 27 Atoms | 90230010 |



MOLECULAR MODEL SET - SODIUM CHLORIDE - OPEN 36 ATOMS

Material : PS/PE

ELEMENTS

SODIUM (VI) - 18
CHLORINE (VI) - 18

LINKS

CONNECTORS - 80

| Product Code | Particulars | HSN Code |
|--------------|--------------------------------|----------|
| 58032 | Sodium Chloride- Open 36 Atoms | 90230010 |



SINCE 1966

MOLECULAR MODEL SET - SODIUM CHLORIDE - COMPACT 27 ATOMS

Material : PS/PE


ELEMENTS

SODIUM (VI) - 13
CHLORINE (VI) - 14

LINKS

CONNECTORS - 60

| Product Code | Particulars | HSN Code |
|--------------|-----------------------------------|----------|
| 58033 | Sodium Chloride- Compact 27 Atoms | 90230010 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

NEW



MOLECULAR MODEL SET - DIAMOND

Material : PS/PE

ELEMENTS


CHLORINE (IV) - 30

LINKS

CONNECTORS - 45

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 58034 | Diamond | 90230010 |

NEW



MOLECULAR MODEL SET - WATER/ICE (H₂O)

Material : PS/PE

ELEMENTS

OXYGEN (IV) - 26

HYDROGEN (II) - 52

LINKS


CONNECTORS - 110

LINK REMOVER TOOL - 1

| Product Code | Particulars | HSN Code |
|--------------|------------------------------|----------|
| 58035 | Water/Ice (H ₂ O) | 90230010 |

SINCE 1966

NEW



MOLECULAR MODEL SET - SODIUM ZINC SULPHIDE

Material : PS/PE

ELEMENTS

METAL (IV)-22

SULPHUR (IV) - 23

LINKS

CONNECTORS - 70

| Product Code | Particulars | HSN Code |
|--------------|---------------|----------|
| 58036 | Zinc Sulphide | 90230010 |



ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

MOLECULAR MODEL SET - SURCOSE



Material : PS/PE


ELEMENTS

- CARBON (IV) - 12
- HYDROGEN (I) - 22
- OXYZEN (II) - 11

LINKS

- CONNECTORS - 70
- LINK REMOVERTOOL - 1

| Product Code | Particulars | HSN Code |
|--------------|--|----------|
| 58037 | Surcose (C ₁₂ H ₂₂ O ₁₁) | 90230010 |



MOLECULAR MODEL SET - GLUCOSE



Material : PS/PE

ELEMENTS

- CARBON (IV) - 12
- HYDROGEN (I) - 24
- OXYZEN (II) - 12

LINKS

- CONNECTORS - 70
- LINK REMOVERTOOL - 1

| Product Code | Particulars | HSN Code |
|--------------|--|----------|
| 58038 | Glucose (C ₆ H ₁₂ O ₆) | 90230010 |



SINCE 1966

MOLECULAR MODEL SET - FAT (GLYCERYL TRISTERATE) (C₅₇H₁₁₀O₆)

Material : PS/PE


ELEMENTS

- CARBON (IV) - 57
- HYDROGEN (I) - 110
- OXYZEN (II) - 6

LINKS

- CONNECTORS - 200
- LINK REMOVERTOOL - 1

| Product Code | Particulars | HSN Code |
|--------------|---|----------|
| 58039 | Fat (Glyceril Tristerate) (C ₅₇ H ₁₁₀ O ₆) | 90230010 |




✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

NEW



MOLECULAR MODEL CHOLESTEROL (C₂₇H₄₆O)

Material : PS/PE

ELEMENTS

- CARBON (IV) - 27
- HYDROGEN (I) - 46
- OXYZEN (II) - 1

LINKS

- CONNECTORS - 110
- LINK REMOVERTOOL - 1

| Product Code | Particulars | HSN Code |
|--------------|---|----------|
| 58041 | Cholesterol (C ₂₇ H ₄₆ O) | 90230010 |

NEW



MOLECULAR MODEL SET - SULPHUR

Material : PS/PE

ELEMENTS

- SULPHUR (II) - 24

LINKS

- CONNECTORS - 45
- LINK REMOVERTOOL - 1

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 58041 | Sulphar | 90230010 |

SINCE 1966

NEW



MOLECULAR MODEL SET - BUCKMINSTER FULLERENE - "BUCKY BALL" - C₆₀

Material : PS/PE

ELEMENTS

- CARBON (IV) - 60

LINKS

- CONNECTORS - 100

| Product Code | Particulars | HSN Code |
|--------------|---|----------|
| 58042 | Buckminster Fullerence "Bucky Ball" - C ₆₀ | 90230010 |

✉ info@asianhyd.com





ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

MOLECULAR MODEL SET - MOLECULAR GEOMETRY SET (VSEPR)

Material : PS/PE

ELEMENTS

- BORON (VI) - 1
- CARBON (II,IV) - 7
- CHLORINE(II,IV) - 12
- FLUORINE (I,IV) - 15
- HYDROGEN (I) - 8
- METAL (II,IV) - 3
- NITROGEN (I) - 6
- OXYGEN (I,IV) - 9
- PHOSPHORUS (V) - 1
- SULPHUR (II,IV,VI) - 4
- XENON (V) - 1

LINKS

- CONNECTORS - 85
- ORBITAL 2D - 15
- LINK REMOVER TOOL - 1



| Product Code | Particulars | HSN Code |
|--------------|--|----------|
| 58043 | Molecular Geometry (VSEPR) - (Valence Shell Electron Pair Repulsion) Set | 90230010 |

MOLECULAR MODEL SET - SHAPES OF MOLECULAR - 8 MODEL KIT

Material : PS/PE


Contains sufficient part to make the eight atomic models shapes (HCl, BeCl₂. H₂O, Bh₃, Nh₃, Ch₄, PCI₅, Sf₆). Each model is an example of different orientation of the bonds and cover coordination numbers 1-6. Lone pairs are represented by brown spheres or brown pear shaped parts.

ELEMENTS

- BORON (III) - 1
- CARBON (IV) - 1
- CHLORINE(I,IV) - 8
- FLUORINE (I) - 9
- HYDROGEN (I) - 13
- METAL (I) - 1
- NITROGEN (IV) - 6
- OXYGEN (IV) - 9
- PHOSPHORUS (V) - 1
- SULPHUR (VI) - 4
- SPHERICAL LONE PAIR -
- ORBITAL(I) - 6

LINKS

- CONNECTORS - 6
- PEAR SHAPED BONDS - 32



| Product Code | Particulars | HSN Code |
|--------------|----------------------------------|----------|
| 58044 | Shapes of Molecules- 8 Model Kit | 90230010 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS



MOLECULAR MODEL SET - MOLECULAR ATOMIC ORBITAL SET - 14 MODEL KIT


Material : PS/PE

This set contains enough individual parts to build 14 distinct models

- 1s Unhybridised (1)
- 2s Unhybridised (1)
- 2s Unhybridised (3)
- 3d Unhybridised (5)
- 2s plus three 2p orbitals Unhybridised (1)
- sp, sp², sp³ Hybridised (3)

Pink and Purple colored pear shaped lobes represent the two wave phase (+ve and -ve) of the p and d atomic orbitals. White sphere represents the atomic nuclei. Each model comes with its own transparent base for display.

| Product Code | Particulars | HSN Code |
|--------------|---------------------------------|----------|
| 58045 | Atomic Orbital Set 14 Model Kit | 90230010 |



CRYSTAL MODEL SET

Material : Polyethylene/LDPE

These Crystal Model Sets contain different coloured jacks (Red, Black & White) & connectors (long & short) in a prescribed number to be joined together in an illustrated way to form structure like Diamond, Graphite & Sodium Chloride.

| Product Code | Particulars | HSN Code |
|--------------|-----------------|----------|
| 58101 | Diamond | 90230010 |
| 58102 | Sodium Chloride | 90230010 |
| 58103 | Graphite | 90230010 |



ASIAN SCIENTIFIC INSTRUMENTS
SUPER MARKET FOR LAB ITEMS

COPLIN JAR (EURO DESIGN)

Material : Polyethylene

These Coplin Jar made of polypropylene can hold 10 microscope slides of the size 75mm x 25mm. Slides are held upright due to the grooved interior. The built in seal ring in the flat cap makes this Coplin Jar a completely leak proof one..



| Product Code | Particulars | HSN Code |
|--------------|--------------------------|----------|
| 61101 | Coplin Jar (Euro Design) | 39231090 |

COPLIN JAR

Material : Polyethylene/

Polylab Coplin Jar made of Polypropylene can hold 10 slides of the size 1"x3" back to back. The interior of the Coplin Jar is grooved in a way to hold slides in erect position. Domed & shallow thread screw cap provides protection to the slides.




| Product Code | Particulars | HSN Code |
|--------------|--------------------------|----------|
| 62101 | Coplin Jar (Euro Design) | 39231090 |

SINCE 1966

SLIDE MAILER

Material : Polyethylene

Slide Mailers present the best option to carry prepared slides from one place to another. Moulded in Polypropylene, these Slide Mailers are available in four sizes i.e. 1, 2, 5 & 10 slides.



| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 62100 | 1 Slide | 39231090 |
| 62102 | 2 Slide | 39231090 |
| 62103 | 5 Slide | 39231090 |
| 62111 | 10 Slide | 39231090 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS




SLIDE MAILER (TRANSPARENT)

Material : PP Material Conforming to US FDA 21CFR.

These Transparent Slide mailers made of polypropylene are the best suited option for secure shipping & storage of slides. Available in different sizes, these slide mailers provide an excellent view of the slides placed insides for easy identification.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 62201 | 1 Slide | 39231090 |
| 62202 | 2 Slide | 39231090 |
| 62203 | 3 Slide | 39231090 |
| 62204 | 5 Slide | 39231090 |



SLIDE BOX

Material : Polystyrene

These Slide Boxes are durable, compact & provide utmost protection to 1"x3" slides. Heavy walls will not warp, splinter or crack. Boxes are unaffected by humidity and are thoroughly insect proof. base of the slide box has a numbered index for easy slide identification.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 62104 | 25 Slide | 39231090 |
| 62105 | 50 Slide | 39231090 |
| 62106 | 100 Slide | 39231090 |



ASIAN SCIENTIFIC INSTRUMENTS
SUPER MARKET FOR LAB ITEMS

SLIDE BOX WITH HINGES

Material : Polystyrene

Adding to the already existing range, these slide boxes are equipped with hinges & lock which makes them the perfect option for keeping the slides safe & secure.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 62211 | 25 Slide | 39231090 |
| 62212 | 50 Slide | 39231090 |
| 62113 | 100 Slide | 39231090 |




SLIDE TRAY

Material : Polystyrene

Moulded in high impact polystyrene this slide tray is an excellent option for holding upto 20 microscopic slides at an angle for easy handling. Placement of slides has been designed in such a way that slides remain protected even when the trays are stacked .

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 62301 | 20 Slide | 39269099 |




SINCE 1966

SLIDE FILE

Material : Polystyrene/SAN

Polylab Slide File can hold upto a maximum of 20 microscopic slides.The base moulded in HIPS horizontally holds the slide.The spaces are numbered for easy identification of the slides.The transparent cover moulded in SAN provides easy reading of the ID labels.


| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 62351 | 20 Slide | 39269099 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS




SLIDE STORAGE RACK

Material : Acrylonitrile Butadiene Styrene/SAN

Adding to the existing range of slide accessories, this slide storage rack made of ABS can hold 100 slides, recognizable through the numeric indexing. Spacing between the slots provide adequate backward & forward movement of slides for convenient recognition. The transparent cover provide excellent view of the slides stored inside.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 62501 | 100 Slide | 39231090 |



SLIDE DRAINING TRAY

Material : Acrylonitrile Butadiene Styrene

This slide draining tray moulded in ABS has numeric indexing for easy identification of the slides. Spacing between the slots provide adequate backward & forward movement of the slides and permits two slides per slot thus providing storage for up to 200 slides.

| Product Code | Particulars | HSN Code |
|--------------|----------------|----------|
| 62502 | 100/200 Slides | 39269099 |

SINCE 1966



FLOAT RACK

Material : Polypropylene

This Float Rack moulded in Polypropylene is quite compact, autoclavable and can hold 16 micro centrifuge tubes for incubation in water baths. Easy to load, these are essential in microbiology & immunology labs. The surface of the rack has raised alpha-numeric index for easy identification and has a central rod handle for convenient lifting.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 63201 | Float Rack | 39269099 |



ASIAN SCIENTIFIC INSTRUMENTS
SUPER MARKET FOR LAB ITEMS

THREE TIER MICRO TUBE RACK

Material : Polypropylene

The interlocking design of these 3 tier micro tube racks makes them ideal for caring a user friendly work station. Unique tired design can carry 3 different size of micro tubes on each tier simultaneously and makes sample loading and transfer an easy process.

Bottom row carries 24x0.2ml (2x12 array) micro tubes

Middle row carries 14x0.5ml (2x7array) micro tubes

Upper row carries 12x1.5ml/2.0ml (2x6 array) micro tubes

| Product Code | Particulars | HSN Code |
|--------------|------------------------|----------|
| 63205 | 3 Tier Micro Tube Rack | 39269099 |



MCT Twin Rack

Material : Polypropylene

PolyLab MCT Twin racks are designed to hold 0.5, 1.5, 2.0 ml micro centrifuge tubes. The transparent top provides an excellent view as well as protection cover for the stored tubes. Each face of the base has a unique alphanumeric index for easy identification of the samples.

| Product Code | Particulars | HSN Code |
|--------------|---|----------|
| 63211 | Rack for 0.5ml & 1.5ml MCTS - 108 place (Holds 108, 1.5/2.0ml tubes where as the reverse side holds 108, 0.5ml tubes. | 39231090 |
| 63215 | Rack for 0.5ml & 1.5ml MCTS - 80 place (Holds 80, 1.5/2.0ml tubes where as the reverse side holds 60, 0.5ml tubes. | 39231090 |




MCT RACK

Material : Polypropylene

These racks holds 80, 1.5ml/2.0ml tubes (5x16array). Embossed alphanumeric index provides quick sample identification and retrieval.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 63216 | 20 Slide | 39231090 |



SINCE 1966



✉ info@asianhyd.com


ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS


PCRTUBE RACK

Material : Polycarbonate/Polypropylene

This rack, as the name suggests, is used for keeping PCR Tubes. Top & base of this rack are made of tough & transparent polycarbonate while the middle portion is made of polypropylene. The tube holding spaces are uniformly placed so as to facilitate the use of multi pipettes. A no. of racks can be stacked together to save shelf space.



| Product Code | Particulars | HSN Code |
|--------------|--------------------------------|----------|
| 63212 | Rack for 96 PCR Tubes of 0.2ml | 39231090 |



PCR TUBE RACK WITH HINGES-AUTOCLAVABLE

Material : Polyethylene

These autoclavable PCR racks can accommodate individual 0.2ml PCR tubes and strips, have hinged lid and alphanumeric identification for quick location and retrieval of tubes. These racks are perfect for pre or post PCR sampling.

| Product Code | Particulars | HSN Code |
|--------------|----------------------------|----------|
| 63217 | PCR Tube Racks with Hinges | 39231090 |

✉ info@asianhyd.com






ASIAN SCIENTIFIC INSTRUMENTS SUPER MARKET FOR LAB ITEMS

MCT BOX

Material : Polycarbonate

MCT Boxes are moulded in polycarbonate and are therefore strong as well as autoclavable. One size can hold 64MCTs of 1.5ml & the other can hold 100MCT's of 0.5ml. These boxes are designed to work at temperatures ranging from -190C to -1210C. The transparent cover provides a good view of the contents inside & is numbered for better inventory control.

| Product Code | Particulars | HSN Code |
|--------------|----------------------------|----------|
| 63213 | Rack for 100 MCTs of 0.5ml | 39231090 |
| 63214 | Rack for 64 MCTs of 1.5ml | 39231090 |



CENTRIFUGE TUBE CONICAL BOTTOM

Material : PP Material Conforming to US FDA 21 CFR.

Polylab Centrifuge tubes with conical bottom are available in two sizes i.e. 15ml & 50ml. Moulded in Polypropylene, these tubes are clear, strong & safe. These autoclavable tubes are provided with printed graduations & writing area along with length on the outer surface. The built in seal ring in the cap makes these tubes absolutely leakproof. These tubes can withstand forces upto 7500xG & 12000xG for 15ml. & 50ml. tubes respectively.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 63301 | 15ml | 39269099 |
| 63302 | 50ml | 39269099 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS
SUPER MARKET FOR LAB ITEMS



RACK FOR CENTRIFUGE TUBE (FOLDING) SPACE SAVER

Material : Polypropylene

These space saver folding centrifuge tube racks are molded in PP. Have raised tube sockets which guarantees firm hold of the tubes at their exact position. The alpha numeric grid reference on the top tier ensures easy tube location & sample positioning. They are supplied assembled and folded, can be ready for use in few seconds with a simple action and can be folded again after use for reducing storage area. These racks are autoclavable & can withstand temperature ranging between -80°C to +121°C. These racks are available in individual for 15ml to 50ml tubes and in combination for both tubes.

| Product Code | Particulars | HSN Code |
|--------------|---|----------|
| 63341 | Space Saver 15ml - 45 Place | 39269099 |
| 63342 | Space Saver 50ml - 18 Place | 39269099 |
| 63343 | Space Saver 15ml & 50ml - 21 Place & 12 Place resp. | 39269099 |

SINCE 1966



CENTRIFUGE TUBE BOX

Material : Polypropylene

Polylab Centrifuge tube boxes are molded in polypropylene. These boxes provide excellent substitute for conventional cardboard or foam racks. The transparent top gives a clear view of the tubes placed in the box which are easily identifiable because of the indexed grid. These boxes are autoclavable and can withstand temperatures ranging from -80°C to +121°C.

| Product Code | Particulars | HSN Code |
|--------------|--------------------|----------|
| 63351 | For 15ml-36 Places | 39231090 |
| 63352 | For 50ml-16 Places | 39231090 |

✉ info@asianhyd.com



ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

CENTRIFUGE TUBE ROUND BOTTOM

Material : Polypropylene/Polycarbonate

PolyLab Centrifuge Tubes with round bottom are moulded in Polypropylene & Polycarbonate these tubes have excellent clarity and are strong enough to be used upto 50,000xG. These tubes can be autoclaved, however, repeated autoclaving in PC tubes reduce their mechanical strength.

| Product Code | Particulars | HSN Code |
|--------------|--------------|----------|
| 63401 | 16ml (P.P) | 39269099 |
| 63402 | 50ml (P.P) | 39269099 |
| 63411 | 16ml (P.C) | 39269099 |
| 63412 | 50ml (P.C) | 39269099 |
| 63421 | Cap for 16ml | 39235010 |
| 63422 | Cap for 50ml | 39235010 |



OAK RIDGE CENTRIFUGE TUBES

Material : Polypropylene

Oak Ridge Centrifuge Tubes, made of Polypropylene, have excellent contact clarity and can be used upto 50,000xG. The screw cap is fitted with silicon O-Ring making it totally leakproof. These tubes should not be used below 4°C.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 63501 | 30ml | 39269099 |
| 63502 | 50ml | 39269099 |
| 63503 | 100ml | 39269099 |

Actual brim capacity of 63502 is 45ml.



SINCE 1966

TEST TUBE WITH SCREW CAP

Material : Polypropylene

These polypropylene tubes provide excellent option for RIA, coagulation & bacteriology etc. The cap is screwed on to the tube & hence this tube can be filled up to the brim

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 64016 | 16x100mm | 39269099 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS



CAP FOR TEST TUBE

Material : Polypropylene

This transparent autoclavable polypropylene cap for Test Tubes proves to be an excellent protection to the contents to the tubes. The grooves inside the cap tighten the fitting on the tube while the outer grooves would hold identification labels. The built in seal ring would set in the tube making the combination a leakproof one.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 64021 | Ø 25mm | 39235010 |



PLASTIC STOPPER

Material : Polypropylene

Polylab Plastic Stoppers made in polypropylene provide excellent option for replacing glass stoppers for bottles & flasks. Top knob has been designed in such a way that handling the stopper is pain free exercise. These stoppers when placed on the bottle or flask neck provide completely leak proof combination.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 64031 | 10/19 | 39235010 |
| 64032 | 12/21 | 39235010 |
| 64033 | 14/23 | 39235010 |
| 64034 | 19/26 | 39235010 |
| 64035 | 24/29 | 39235010 |
| 64036 | 29/32 | 39235010 |

SINCE 1966



STORAGE VIAL

Material : PP Material Conforming to US FDA 21 CFR.

As the name suggests, these vials are commonly used for storage of biological material, human & animal cells. The Cap is provided with a built in seal ring making the vial leakproof. The base of the vial has longitudinal grooves which give it a non-twisting character when placed in a vial rack.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 64102 | 2.0ml | 39233090 |
| 64102 | 5.0ml | 39233090 |



ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

COLLECTION VIAL

Material : PP Material Conforming to US FDA 21CFR.

Moulded in PP these autoclavable collection vials can be used for storage of different PCT reagents, enzymes & other diagnostics & biochemical reagents & samples. The inbuilt seal ring in the cap makes these vials leak proof.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 64301 | 5ml | 39233090 |



SCINTILLATION VIAL

Material : High Density Polyethylene

These unbreakable vials, moulded in Polyethylene, are fitted with leakproof screw caps. These vials can be used with standard counting equipment for liquid scintillation system and are excellent for low activity count.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 65001 | 8ml | 39233090 |
| 65002 | 20ml | 39233090 |



SINCE 1966

RACK FOR SCINTILLATION VIAL

Material : Polypropylene

These Racks for Scintillation Vials, moulded in Polypropylene, are autoclavable & corrosion resistant. There is an individual rack for both 8ml. & 20ml. vials. These rugged & stable racks can hold scintillation vials for convenient storage, handling & sample preparations. Can be stacked when not filled for easy storage. Embossed alphanumeric index provides quick sample identification & retrieval.

| Product Code | Particulars | HSN Code |
|--------------|------------------|----------|
| 65101 | 90 vials of 8ml | 39269099 |
| 65101 | 50 vials of 20ml | 39269099 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

CRYO VIAL

Material : PP Material Conforming to US FDA 21CFR.

Cryo vials are normally used for storage of Biological material, human & animal cells. These sterilized tubes are provided with white , prominent graduations and a long marking area for sample identification. The cap of these vials is provided with a seal ring inside to make it completely leakproof. The longitudinal grooves in the base give a non-twisting features to these vials when placed in a cryo rack.



| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 66002 | 1.8ml | 39233090 |
| 66003 | 4.5ml | 39233090 |



CRYO RACK

Material : Polycarbonate

Polylab Cryo Rack is moulded in Polycarbonate making it tough, durable & autoclavable. This handy rack can hold 50 vials at a time. The unique feature of this rack is that vials can be unscrewed without actually picking them up as the vial will lock in every well and hence will not twist. This rack is provided with four anti skid rubbers & each well is identified with an alpha-numeric index.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 66201 | Cryo Rack | 39233090 |

✉ info@asianhyd.com





ASIAN SCIENTIFIC INSTRUMENTS SUPER MARKET FOR LAB ITEMS

CRYO BOX (P.C)

Material : Polycarbonate

Polylab Cryo Boxes are moulded in Polycarbonate making it extra strong and autoclavable. These can accommodate Cryo Vials of up to 4.5ml. capacity and are designed to be used at temperatures ranging from -196°C + 121°C. The transparent cover of the box gives a good view of the contents inside and is numbered for better inventory control. The specially designed venting in base & cover ensures a smooth flow for vapours for consistency in cooling. Dimensions of these boxes are designed to be compatible with vertical, horizontal & drawers style freezer racks.

| Product Code | Particulars | HSN Code |
|--------------|-----------------------------|----------|
| 66300 | 50 Places for 1ml or 1.8ml | 39231090 |
| 66301 | 81 Places for 1ml or 1.8ml | 39231090 |
| 66302 | 81 Places for 4.5ml | 39231090 |
| 66303 | 100 Places for 1ml or 1.8ml | 39231090 |



CRYO BOX (P.P)

Material : PP Material Conforming to US FDA 21CFR.

Moulded in Polypropylene, the top of the box is attached to the bottom through well designed built-in hinges and the locking systems ensures secure closer preventing accidental opening of these boxes in deep freezers. They are stackable & designed for sample storage up to -80°C. The bottom has built-in spaces to hold vials up to 2ml. capacity. Available in two different sizes, these autoclavable boxes bear alpha numeric index for easy identification of different vials placed inside.


| Product Code | Particulars | HSN Code |
|--------------|-----------------------------|----------|
| 66451 | 50 Places for 1ml or 1.8ml | 39231090 |
| 66501 | 81 Places for 1ml or 1.8ml | 39231090 |
| 66502 | 100 Places for 1ml or 1.8ml | 39231090 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS



MICRO TIP BOX

Material : PP Material Conforming to US FDA 21CFR.

Polylab Micro tip boxes molded in high quality polypropylene can be autoclavable. The unique feature of these tip boxes is the robust hinge design which prevents the top from being easily disengaged from the bottom of the box and consistent wall thickness for repetitive use in single and multi channel pipette. Transparent top cover provides excellent view of the micro tips inside the box.

| Product Code | Particulars | HSN Code |
|--------------|----------------------------------|----------|
| 67011 | For 96 Tips of 2-10 μ l | 39231090 |
| 67012 | For 96 Tips of 2-200 μ l | 39231090 |
| 67013 | For 100 Tips of 200-1000 μ l | 39231090 |



MICRO TIP BOX (EURO DESIGN)

Material : PP Material Conforming to US FDA 21CFR.

Adding to the previous range of Polylab Micro tip boxes, these boxes are designed with uniform tip holding space to facilitate easy use in multichannel pipette. All tip boxes contains 96 tip holes.

| Product Code | Particulars | HSN Code |
|--------------|---------------------------------|----------|
| 67121 | For 96 Tips of 2-10 μ l | 39231090 |
| 67122 | For 96 Tips of 2-200 μ l | 39231090 |
| 67123 | For 96 Tips of 200-1000 μ l | 39231090 |

SINCE 1966

VOLUMETRIC FLASK

Material : Polypropylene

Polylab Volumetric Flasks, Blown in polypropylene are strong, unbreakable, individually calibrated and have good contact clarity. These autoclavable flasks are provided with screw caps which makes them completely leak proof.



| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 68001 | 25ml | 39233090 |
| 68002 | 50ml | 39233090 |
| 68003 | 100ml | 39233090 |
| 68004 | 250ml | 39233090 |
| 68005 | 500ml | 39233090 |
| 68006 | 1000ml | 39233090 |



ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

PETRI DISHES DISPOSABLE GAMMA STERILE

Material : PS

Polylab disposable petri dishes are made of PS confirming to US FDA 21 CFR. Its dimensional accuracy of 90mm ensures compatibility with most automated systems. Its optical clarity ensures a clear view on microscope. These triple vented dishes have elevated lid that allows for good gas exchange & encourage evaporation. These petri dishes are sterilized by radiation.

| Product Code | Particulars | HSN Code |
|--------------|--|----------|
| 69001-G | 90mm(Individually Packed) - Gamma Sterile | 39269099 |
| 69002-G | 90mm(Sleeve of 20) - Gamma Sterile | 39269099 |



PETRI DISHES DISPOSABLE ETO STERILE

Material : PS

Apart from all the feature of disposable petri dishes are sterilized by ETO.

| Product Code | Particulars | HSN Code |
|--------------|--|----------|
| 69001-E | 90mm(Individually Packed) - ETO Sterile | 39269099 |
| 69002-G | 90mm(Sleeve of 20) - ETO Sterile | 39269099 |



PETRI DISHES DISPOSABLE ASEPTIC

Material : PS

These petri dishes are manufactured in clean room under aseptic conditions which ensures no contact with any contaminants and best performance in all application. It also eliminates discoloration & stress caused by radiation sterilization.

| Product Code | Particulars | HSN Code |
|--------------|--|----------|
| 69001-A | 90mm(Individually Packed) - Aseptic | 39269099 |
| 69002-A | 90mm(Sleeve of 20) - Aseptic | 39269099 |



RACK FOR PETRI DISHES

Material : Polymethyl Methacrylate/ Polycarbonate

These Racks are clear acrylic construction with polycarbonate posts. 69101 has a capacity to hold 60 petri dishes of 90mm size while 69102 one can hold 56 petri dishes of 60mm. These racks are very useful during inoculation, incubation & storage. Visibility from all sides enables the culture to be seen & checked during incubation.

| Product Code | Particulars | HSN Code |
|--------------|---------------|----------|
| 69101 | Rack for 90mm | 39269099 |
| 69102 | Rack for 60mm | 39269099 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS




UNIVERSAL MULTI RACK

Material : Polypropylene

Universal multi rack is an excellent lab apparatus compatible with different sizes of tubes. This autoclavable rack may be rotated to any of the four positions.

- 1st side - 8 wells with \varnothing 30mm for 50ml tubes
- 2nd side - 24 wells with \varnothing 17.5mm for 15ml tubes
- 3rd side - 20 wells with \varnothing 20mm
- 4th side - 40 wells with \varnothing 12mm for RIAVIAL

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 76001 | Multi Rack | 39269099 |



UNIVERSAL COMBI MULTI RACK

Material : Polypropylene

This interlocking racks are used as stand alone or infinite units can be linked to customized your specific requirements. Four sides of these autoclavable racks can accommodate different sizes of tubes.

- 1st side - 4 wells with \varnothing 30mm for 50ml tubes
- 2nd side - 32 wells with \varnothing 8.2mm for 0.5ml tubes
- 3rd side - 15 wells with \varnothing 17.5mm
- 4th side - 32 wells with \varnothing 11.6mm for 1.5/2.0ml tubes

| Product Code | Particulars | HSN Code |
|--------------|-----------------|----------|
| 76003 | Combi Link Rack | 39269099 |

SINCE 1966



ROTATABLE MULTI RACK

Material : Polypropylene

This autoclavable rotatable multi rack moulded in polypropylene is a combination of four modules providing spaces for different sizes of test tubes and centrifuge tubes. These modules can be rotated in different directions to accommodate different types and sizes of tubes at the same time.

- 1st side - 2 wells with \varnothing 30mm for 50ml tubes
- 2nd side - 10 wells with \varnothing 13.5mm for 10-13ml tubes
- 3rd side - 8 wells with \varnothing 20mm for 10-15ml tubes
- 4th side - 6 wells with \varnothing 11.6mm for 15-20ml tubes

| Product Code | Particulars | HSN Code |
|--------------|----------------------|----------|
| 76002 | Rotatable Multi Rack | 39269099 |

✉ info@asianhyd.com





ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

ROTATABLE MULTI RACK-EXPANDABLE

Material : Polypropylene

Adding to the existing rotatable multi rack these racks has a unique rotate and lock design that helps add or subtract racks to create configurations to fit number of wells and tubes according to user's needs.


Supplied in set of 4 units (in assorted colors). Each unit has 4 sides with different hole sizes.

1st side accommodates 2 units of 50ml centrifuge tubes

2nd side accommodates 8 units of 15ml centrifuge tubes

3rd side accommodates 9 units of 1.5/0.2ml MCT and 6 units of 0.5ml MCT.

4th side accommodates 48 units of 0.2ml PCR or 6 units of 8 tube PCR strips or 8 units of 12 tubes PCR strips (when 2 units locked together) or 1 unit of 96 well (unskirted or semi skirted) PCR plate (white 2 units locked together)




| Product Code | Particulars | HSN Code |
|--------------|---------------------------------|----------|
| 76004 | Rotatable Multi Rack-Expandable | 39269099 |

NESTLER CYLINDER STAND

Material : Polypropylene

This stand, as the name suggests, is used for keeping nestler cylinders. Made of polypropylene, this stand is available in two sizes, one for 50ml nestler cylinder & other for 100ml, nestler cylinder.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 76101 | 50ml | 39269099 |
| 76102 | 100ml | 39269099 |



SINCE 1966


TEST TUBE STAND

Material : Polypropylene

This autoclavable stand is perhaps the most commonly used Test tube Stand in various laboratories. For better visibility, this 6 place rack is end supported by 2 columns. The top plate has 2 holes for holding 25mm dia test tubes, four holes to hold 16mm dia tubes where as the base has corresponding hemispherical wells. It also has six vertical pins for drying tubes.

Another option available in this section is a test tube stand in which all the holes in the top plates can hold 25mm dia tubes.

| Product Code | Particulars | HSN Code |
|--------------|----------------------------|----------|
| 77701 | 16mm & 25mm θ Tubes | 39269099 |
| 77702 | 25mm θ Tubes | 39269099 |



✉ info@asianhyd.com


ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

TEST TUBE STAND (ROUND)

Material : Polypropylene

A handy, little, space saver, this circular twelve place Polypropylene Test Tube Stand has four 25mm dia holes & eight 19mm dia holes in its top plate. The base has matching hemispherical wells & 16 vertical pins for drying tubes surrounded by trough to catch drippings. This autoclavable rack can be easily disassembled for cleaning for cleaning purpose.



| Product Code | Particulars | HSN Code |
|--------------|----------------------------|----------|
| 77704 | 19mm & 25mm θ Tubes | 39269099 |



RACK FOR MICRO CENTRIFUGE TUBES

Material : Polypropylene

Moulded in Polycarbonate, these autoclavable racks are made for holding 1.5ml. micro centrifuge tubes. These racks have a unique stacking feature which prevents movement in any direction. The holes are numbered & lettered for easy identification.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 77705 | 24 Tubes | 39269099 |
| 77705A | 48 Tubes | 39269099 |

✉ info@asianhyd.com





ASIAN SCIENTIFIC INSTRUMENTS


SUPER MARKET FOR LAB ITEMS

TEST TUBE STAND (3 TIER) P.P

Material : RPP

These ingenious racks with a 3-tier design provide clear view of the tube contents & convenient insertion & removal of tubes. These autoclavable Racks can withstand sub-freezing temperatures.

| Product Code | Particulars | HSN Code |
|--------------|-----------------|----------|
| 77706 | 13mm x 31 Tubes | 39269099 |
| 77707 | 16mm x 31 Tubes | 39269099 |
| 77708 | 20mm x 20 Tubes | 39269099 |
| 77709 | 25mm x 12 Tubes | 39269099 |
| 77710 | 25mm x 18 Tubes | 39269099 |
| 77711 | 32mm x 12 Tubes | 39269099 |
| 77712 | 25mm x 36 Tubes | 39269099 |
| 77714 | 13mm x 18 Tubes | 39269099 |
| 77715 | 15mm x 18 Tubes | 39269099 |
| 77716 | 18mm x 12 Tubes | 39269099 |
| 77717 | 13mm x 62 Tubes | 39269099 |
| 77718 | 16mm x 62 Tubes | 39269099 |
| 77719 | 20mm x 40 Tubes | 39269099 |
| 77720 | 13mm x 48 Tubes | 39269099 |
| 77721 | 13mm x 72 Tubes | 39269099 |
| 77722 | 25mm x 24 Tubes | 39269099 |




TEST TUBE STAND (3 TIER) P.C

Material : Polycarbonate

These polycarbonate test tube stands features ergonomic designed for easy handling. Provides good chemical resistance. These autoclavable stands can with stands temperature from -135oC to +135oC.

| Product Code | Particulars | HSN Code |
|--------------|-----------------|----------|
| 77801 | 13mm x 31 Tubes | 39269099 |
| 77802 | 16mm x 31 Tubes | 39269099 |
| 77803 | 20mm x 20 Tubes | 39269099 |
| 77804 | 25mm x 12 Tubes | 39269099 |
| 77805 | 13mm x 62 Tubes | 39269099 |
| 77806 | 16mm x 62 Tubes | 39269099 |
| 77807 | 20mm x 40 Tubes | 39269099 |
| 77808 | 25mm x 24 Tubes | 39269099 |
| 77809 | 25mm x 18 Tubes | 39269099 |
| 77810 | 25mm x 36 Tubes | 39269099 |
| 77811 | 32mm x 12 Tubes | 39269099 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS




TEST TUBE RACK

Material : Polypropylene

These polypropylene assembled test tube racks are an excellent option for storing large number of samples. Manufactured from robust, sturdy and lightweight material makes it suitable for easy movement and transport. These autoclavable rack are available in three sizes to fit standard 10-13mm or 13-16mm or 18-20mm tubes.

| Product Code | Particulars | HSN Code |
|--------------|-----------------|----------|
| 77706 | 13mm x 31 Tubes | 39269099 |
| 77707 | 16mm x 31 Tubes | 39269099 |
| 77708 | 20mm x 20 Tubes | 39269099 |
| 77709 | 25mm x 12 Tubes | 39269099 |
| 77710 | 25mm x 18 Tubes | 39269099 |
| 77711 | 32mm x 12 Tubes | 39269099 |
| 77712 | 25mm x 36 Tubes | 39269099 |
| 77714 | 13mm x 18 Tubes | 39269099 |
| 77715 | 15mm x 18 Tubes | 39269099 |
| 77716 | 18mm x 12 Tubes | 39269099 |
| 77717 | 13mm x 62 Tubes | 39269099 |
| 77718 | 16mm x 62 Tubes | 39269099 |
| 77719 | 20mm x 40 Tubes | 39269099 |
| 77720 | 13mm x 48 Tubes | 39269099 |
| 77721 | 13mm x 72 Tubes | 39269099 |
| 77722 | 25mm x 24 Tubes | 39269099 |



TEST TUBE PEG RACK

Material : Polycarbonate

These racks can conveniently hold inverted tubes for drying purpose and minimizes collection of airborne contaminants inside the tube. This rack can also be used to hold and dry electrophoresis and chromatography plates.

| Product Code | Particulars | HSN Code |
|--------------|--------------------------|----------|
| 77851 | 96 Places for 13mm Tubes | 39269099 |
| 77852 | 50 Places for 16mm Tubes | 39269099 |




ASIAN SCIENTIFIC INSTRUMENTS SUPER MARKET FOR LAB ITEMS

TEST TUBE STAND (WIRE PATTERN)

Material : Polypropylene

An economical substitute for wire racks, here are these submersible, autoclavable racks which are stackable when empty. These are available in different colors and are made of a special blend of Polypropylene that makes it sink in a water bath and maintain their stability even in agitated water. These space saving racks can easily be disassembled for cleaning purposes. The base is in the form of a grid with square openings and has fold up sides and end bars that lock into rack tops.



| Product Code | Particulars | HSN Code |
|--------------|-----------------|----------|
| 77901 | 13mm x 90 Tubes | 39269099 |
| 77902 | 16mm x 60 Tubes | 39269099 |
| 77903 | 20mm x 40 Tubes | 39269099 |
| 77904 | 25mm x 24 Tubes | 39269099 |
| 77905 | 30mm x 21 Tubes | 39269099 |

TEST TUBE STAND (WIRE PATTERN-FIX)

Material : Polypropylene

These are one piece molded racks and have a three tier grid design. They have good chemical & thermal resistance.



| Product Code | Particulars | HSN Code |
|--------------|-----------------|----------|
| 77906 | 13mm x 36 Tubes | 39269099 |
| 77907 | 16mm x 36 Tubes | 39269099 |

SINCE 1966

RACK FOR MICRO CENTRIFUGE (FOLDING)

Material : Polypropylene

These two-tier versatile racks, made of propylene, are autoclavable & can be used dry, in water baths or in freezers. These racks are provided with alpha numeric index for easy sample identification. Racks are stack able even when loaded and have purpose designed label spaces on each side.



| Product Code | Particulars | HSN Code |
|--------------|---------------------------|----------|
| 77910 | 100 Places for 1.5ml MCTs | 39269099 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS




TEST TUBE BASKETS

Material : Polypropylene

Polylab Test Tube Baskets are made of Polypropylene. These autoclavable, unbreakable and non-corrosive Baskets are ideal for rinsing & draining glassware. These are also useful for storage of test tubes or culture tubes during sterilization or incubation.

| Product Code | Particulars | HSN Code |
|--------------|---------------|----------|
| 78101 | 140x120x110mm | 39269099 |
| 78102 | 160x160x160mm | 39269099 |



TEST TUBE BASKETS (WITH COVER)

Material : Polypropylene

These baskets add to the existing range of test tube baskets. The additional features of these baskets in the lid that has been provided to keep the apparatus placed in it more secure. Molded in PP these baskets are autoclavable, unbreakable & non corrosive

| Product Code | Particulars | HSN Code |
|--------------|---------------|----------|
| 78111 | 110x120x150mm | 39269099 |
| 78112 | 180x170x160mm | 39269099 |
| 78113 | 230x230x230mm | 39269099 |



DRAINING BASKETS

Material : Polypropylene

As the name suggests, this basket with wide net like structure on all sides prove really useful for draining different laboratory apparatus after washing them. Made of polypropylene, this basket can be repeatedly.

| Product Code | Particulars | HSN Code |
|--------------|---------------|----------|
| 78103 | 400x400x100mm | 39269099 |



DRAINING TRAY

Material : Polypropylene

These autoclavable perforated trays can be placed in sink for washing, cleaning and draining of various laboratory glass and plastic apparatus.

| Product Code | Particulars | HSN Code |
|--------------|---------------|----------|
| 78104 | 400x300x100mm | 39269099 |



ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

PIPETTE STAND (HORIZONTAL)

Material : Polypropylene

This detachable Stand can hold 12 pipettes, 6 on each side, in a horizontal position. Moulded in Polypropylene, this autoclavable stand has two sides which are attached to each other with three metallic rods coated with Polypropylene. This stand can be easily disassembled for space saving storage.

| Product Code | Particulars | HSN Code |
|--------------|-----------------|----------|
| 79101 | For 12 Pipettes | 39269099 |



PIPETTE STAND (VERTICAL)

Material : Polypropylene

This is a stable, non-corrosive, detachable stand which can hold 28 pipettes with maximum dia 14mm. The heavy base measures 225mm in diameter and a Polypropylene coated square rod measuring 375 mm in length is screwed into the base. The top disc that can be fixed on the rod has holes for holding pipettes while the lower one has pilot holes where the pipette tips rest and liquids drain down to the base. These discs can be adjusted on the rod according to the length of the pipette. The arrows provided on the discs should be kept parallel to ensure that corresponding holes are being centered.

| Product Code | Particulars | HSN Code |
|--------------|-----------------|----------|
| 79102 | For 28 Pipettes | 39269099 |



SINCE 1966

PIPETTE STAND (94 PIPETTES-ROTARY)

Material : Polypropylene

Moulded in Polypropylene, this unique Stand holds a maximum of 94 pipettes and rotates 360° on a central vertical axis for convenient selection of any particular pipette. The lower rotating disc has four concentric troughs with tapered sides and open bottoms which cradle the pipette tips protectively and permit water drops to run off. The upper disc has a grid pattern concentric with the bottom disc and keeps the pipettes vertical. Liquids draining off the pipettes are retained on the top of the stationary base by the outer ring. The stand can be easily disassembled & reassembled.

| Product Code | Particulars | HSN Code |
|--------------|-----------------|----------|
| 79103 | For 94 Pipettes | 39269099 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

Asian Estate : 4-I-10, Tilak Road Opposite Posnett Bhavan Hyderabad, Telangana-500 001

Ph. : 93930-01939, 90522-23905, 98480-31317, 93472-45096

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS



PIPETTE BOX

Material : Polypropylene

Made of polypropylene, this Pipette Box is designed in such a way that it can hold Pipettes of different lengths. The top of the box screws on the bottom and the height can be adjusted by screwing or unscrewing the top.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 79111 | Pipette Box | 39269099 |



PIPETTE TRAY WITH 4 COMPARTMENTS

Material : Polypropylene

These trays have 4 compartments to hold about 30 pipettes of 1-2-5-10ml capacities, can also be used in drawers as well as on bench top.

| Product Code | Particulars | HSN Code |
|--------------|--------------|----------|
| 79115 | 420x300x30mm | 39269099 |

SINCE 1966



PIPETTE TRAY

Material : PS

As the name suggests, this basket with wide net like structure on all sides prove really useful for draining different laboratory apparatus after washing them. Made of polypropylene, this basket can be repeatedly.

| Product Code | Particulars | HSN Code |
|--------------|---------------|----------|
| 79116 | 400x400x100mm | 39269099 |



ASIAN SCIENTIFIC INSTRUMENTS


SUPER MARKET FOR LAB ITEMS

REAGENT RESERVOIR

Material : Polypropylene

Moulded in Polypropylene, this autoclavable Reagent Reservoir is helpful for multi-channel pipettes to be used into microtiter plates.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 79121 | 75ml | 39269099 |




UNIVERSAL REAGENT RESERVOIR

Material : Polypropylene

This Universal Reagent Reservoir is a graduated single reagent reservoir with maximum volume of 50ml, reservoir with 8 separate numbered 5ml. reagent channels & a reservoir with 12 separate numbered 5ml. reagent channels. made of polypropylene these reagent reservoirs are repeatedly autoclavable and can be used with all multi-channels pipettes.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 79131 | 50ml | 39269099 |



SINCE 1966

STAINING BOX

Material : Polypropylene

Moulded in Polypropylene, this staining Box helps in staining, fixing, de-staining and handling fragile electrophoresis gels and membranes. Box as well as cover will resist temperatures from -700C to +700C. The LPDE plug provides leakproof seal and allows easy liquid drainage.

| Product Code | Particulars | HSN Code |
|--------------|--------------|----------|
| 79301 | 125x125x50mm | 39269099 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS



MEASURING CYLINDERS HEXAGONAL

Material : PP Material Conforming to US FDA 21CFR.

Polylab measuring cylinders moulded in polypropylene are single piece cylinders. These autoclavable cylinders have raised prominent graduations for easy reading & have excellent contact clarity. The hexagonal base provides excellent stability to the cylinder where as the tapered spout works really well while pouring any reagent.

| Product Code | Particulars | HSN Code |
|--------------|-----------------|----------|
| 80031 | 13mm x 31 Tubes | 39269099 |
| 80032 | 16mm x 31 Tubes | 39269099 |
| 80033 | 20mm x 20 Tubes | 39269099 |
| 80034 | 25mm x 12 Tubes | 39269099 |
| 80035 | 13mm x 62 Tubes | 39269099 |
| 80036 | 16mm x 62 Tubes | 39269099 |
| 80037 | 20mm x 40 Tubes | 39269099 |
| 80038 | 25mm x 24 Tubes | 39269099 |



MEASURING JUGS

Material : PP Material Conforming to US FDA 21CFR.

These handy Jugs, moulded in Polypropylene, are clear, autoclavable and have good chemical resistance. These jugs are provided with a handle for easy & convenient lifting and have raised graduations for easy reading.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 81101 | 250ml | 39269099 |
| 81102 | 500ml | 39269099 |
| 81103 | 1000ml | 39269099 |
| 81104 | 2000ml | 39269099 |

SINCE 1966



MEASURING JUGS (EURO DESIGN)

Material : PP Material Conforming to US FDA 21CFR.

Adding to the existing range of Polylab Jugs, these short form transparent jugs, made of polypropylene are autoclavable and have non-drip spout. These jugs have raised graduations, provide excellent contact clarity & are provided with excellent thumb grip on handle.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 81120 | 250ml | 39269099 |
| 81121 | 250ml | 39269099 |
| 81122 | 250ml | 39269099 |
| 81123 | 250ml | 39269099 |
| 81124 | 250ml | 39269099 |
| 81125 | 250ml | 39269099 |



ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

MEASURING JUGS (PRINTED GRADUATION)

Material : PP Material Conforming to US FDA 21CFR.

Like other PolyLab Measuring Jugs, these jugs, also, moulded in polypropylene have excellent clarity as well as very good chemical resistance. The prominent features of these jugs is the easy to read printed graduations. The rectangular handle provides significant advantage for easy pouring and handling. These jugs are autoclavable.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 81131 | 100ml | 39269099 |
| 81132 | 250ml | 39269099 |
| 81133 | 300ml | 39269099 |
| 81134 | 500ml | 39269099 |
| 81135 | 1000ml | 39269099 |
| 81136 | 5000ml | 39269099 |



CONICAL MEASURES

Material : PP Material Conforming to US FDA 21CFR.

These Conical Measures are very good for measuring small quantities of reagents. Moulded in Polypropylene, these autoclavable Measures have a hexagonal base which provides it excellent stability. The Measures are clear, have good chemical resistance and have raised graduation for easy reading.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 81150 | 12ml | 39269099 |
| 81151 | 25ml | 39269099 |
| 81152 | 50ml | 39269099 |
| 81133 | 125ml | 39269099 |
| 81134 | 200ml | 39269099 |



SINCE 1966

SEPARATORY FUNNEL

Material : Polypropylene/PTFE

PolyLab Separatory Funnel is made of polypropylene, giving it excellent contact clarity and strength. This pear shaped, autoclavable separating funnel has a leakproof threaded screw cap. The Polypropylene stopcock, fitted with a self lubricating PTFE plug, is push fitted to the stem of the separatory funnel.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 81201 | 100ml | 39269099 |
| 81202 | 250ml | 39269099 |
| 81203 | 500ml | 39269099 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS



SCOOP

Material : PP Material Conforming to US FDA 21CFR.

Polylab scoops are excellent replacement for metallic ones. Moulded in polypropylene, these scoops are autoclavable and completely eliminate the danger of rusting or corrosion.

| Product Code | Particulars | HSN Code |
|--------------|---------------|----------|
| 81251 | 2ml | 39269099 |
| 81252 | 5ml | 39269099 |
| 81253 | 10ml | 39269099 |
| 81254 | 25ml | 39269099 |
| 81255 | 50ml | 39269099 |
| 81256 | 100ml | 39269099 |
| 81257 | 250ml | 39269099 |
| 81258 | 500ml | 39269099 |
| 81259 | 1000ml | 39269099 |
| 81260 | Scoop - Set/9 | 39269099 |



LONG HANDLE SCOOP

Material : PP Material Conforming to US FDA 21CFR.

PolyLab long handle scoops are of perfect shape and size for taking samples from wide mouth containers. The longer handle of 60ml and 125ml scoops allow the user to scoop the sample from the distance. They can be used as weighing boats as they remain upright when placed on a flat surface.

| Product Code | Particulars | HSN Code |
|---------------|-------------|----------|
| 81271 | 2ml | 39269099 |
| 81272 | 5ml | 39269099 |
| 81273 | 10ml | 39269099 |
| Material : PS | | |
| 81281 | 50ml | 39269099 |
| 81282 | 100ml | 39269099 |
| 81283 | 250ml | 39269099 |



SPATULA

Material : Polypropylene

Moulded in Polypropylene, Spatula consists of a spoon on one side while the other end is a flat construction with a knife edge. This edge is helpful in breaking up crystals of different salts while the spoon is used to draw salts from the containers

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 81301 | 150mm | 39269099 |
| 81302 | 200mm | 39269099 |

✉ info@asianhyd.com





ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

FORCEPS

Material : Polypropylene/HIPS

These polystyrene tweezers forceps are 4" in length, features serrated tips and pointy tapered ends ensures firm grip and its anti slip thumb and finger ridges added control. Ideal for grasping items from surfaces and also used for medical and first aid purpose.

Whereas Polypropylene Scissor-Type Forceps are 5" in length, are chemically compatible with most lab reagents. Double-action hinges provide extra leverage and long life. The jaws have interlocking teeth for holding objects securely. Scissor type handles with catch for locking the forceps. Can be used as clips for chromatography columns. Also serve as excellent tube clamps.



| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 81321 | Tweezer | 39269099 |
| 81322 | Scissor | 39269099 |

SAFETY GOGGLES WITH VENT

Material : PVC/PC/PP

Frame of PolyLab safety goggles is made up of flexible Phthalates free PVC and the 100% optically clear lens is made up of polycarbonate which makes the goggles light weight and more comfortable to wear. It comes with a readjust able band so one size fits all perfectly offering protection against fine dust, droplets, liquid splashes and impact.

It has 4 indirect vents to prevent fogging .It can be comfortably worn over cotact lenses or regular eye glasses. Our safety goggles are ANSI/ISEA Z87.1-2020 certified and comply with standard EN166:201-PErsonal Eye Protection.

These goggles are fully autoclavable.



| Product Code | Particulars | HSN Code |
|--------------|--------------------------|----------|
| 81352 | Safety Goggles with Rent | 90049090 |

STIRRER

Material : Polyprpylene

Made of polypropylene, these Stirrers are very useful for stirring different solutions in the laboratory. These Stirrers prove to be tough and long lasting.

| Product Code | Particulars | HSN Code |
|--------------|------------------|----------|
| 81401 | Ø-10mm x h-250mm | 39269099 |
| 81404 | Ø-6mm x h-150mm | 39269099 |
| 81405 | Ø-7mm x h-200mm | 39269099 |
| 81406 | Ø-7mm x h-250mm | 39269099 |
| 81407 | Ø-7mm x h-300mm | 39269099 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS




POLICEMEN STIRRING ROD

Material : Polypropylene

This, infact, is a multipurpose apparatus which works as a stirrer as well spatula. The fattened ends of this rod form two small spatulas.

| Product Code | Particulars | HSN Code |
|--------------|-----------------|----------|
| 81411 | Ø-6mm x h-245mm | 39269099 |



SOFT LOOPES

Material : Polypropylene

This non-toxic disposable loop is accurate to semi quantitative stand for handling samples.

| Product Code | Particulars | HSN Code |
|--------------|--------------------|----------|
| 81421 | Soft Loopes -1µl | 39269099 |
| 81422 | Soft Loopes -10µl | 39269099 |
| 81423 | Inoculation Needle | 39269099 |



L SHAPED SPREADER

Material : Polypropylene

This sterile, smooth spreader enables even spreading of liquid samples across the surface of a jar plate without grouping & cutting the medium.

| Product Code | Particulars | HSN Code |
|--------------|-------------------|----------|
| 81431 | L Shaped Spreader | 39269099 |



MICRO PESTLE

Material : Polypropylene

This double ended pestle's one end fits 1.5ml. Micro centrifuge tube & the other one fits 0.5ml. Micro centrifuge tube.

| Product Code | Particulars | HSN Code |
|--------------|-----------------------|----------|
| 81441 | for 0.5 and 1.5ml MCT | 39269099 |
| 81442 | for 1.5ml MCT | 39269099 |



ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

PNEUMATIC TROUGH

Material : Polypropylene

Three trough, made of polypropylene, are autoclavable and are available in different sizes. These can be used for a wide range of purpose including storage sterilizing & drying labware etc.

| Product Code | Particulars | HSN Code |
|--------------|-----------------|----------|
| 81501 | Ø-180 x h-60mm | 39269099 |
| 81502 | Ø-200 x h-100mm | 39269099 |
| 81503 | Ø-250 x h-100mm | 39269099 |



PLANTATION POTS

Material : Polypropylene

Polylab Plantation Pots, moulded in polypropylene, are available in four different sizes. These pots have a downward extended rim which provide convenient lifting even if the pot is filled. The bottom of the pot is provided with holes for drain purpose.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 81601 | 100mm | 39269099 |
| 81602 | 150mm | 39269099 |
| 81603 | 200mm | 39269099 |
| 81604 | 250mm | 39269099 |



SINCE 1966

LABORATORY TRAY

Material : Polypropylene

Needless to say, PolyLab Laboratory Tray moulded in polypropylene is a multipurpose tray that can be used for sterilizing, drying glassware, porcelainware & other laboratory instruments. These steam autoclavable trays have tapered walls which prove good for nesting.

| Product Code | Particulars | HSN Code |
|--------------|--------------|----------|
| 81701 | 450x350x75mm | 39269099 |
| 81702 | 375x300x75mm | 39269099 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

Asian Estate : 4-I-10, Tilak Road Opposite Posnett Bhavan Hyderabad, Telangana-500 001

Ph. : 93930-01939, 90522-23905, 98480-31317, 93472-45096

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS



UTILITY TRAY

Material : Polypropylene

Adding to the already existing range of autoclavable laboratories trays, is the Utility Tray and as the name suggests this Polypropylene Tray can be utilized for purposes ranging from drying glassware to sterilizing and to store any laboratory apparatus.

| Product Code | Particulars | HSN Code |
|--------------|---------------|----------|
| 81721 | 375x350x130mm | 39269099 |
| 81722 | 540x435x130mm | 39269099 |




CARRIER TRAY

Material : Polypropylene

This autoclavable Polypropylene Tray has a unique feature i.e. its handle which proves really helpful in carrying reagents or laboratory apparatus from one place to another, single handedly.

| Product Code | Particulars | HSN Code |
|--------------|---------------|----------|
| 81731 | 380x240x115mm | 39269099 |

SINCE 1966



INSTRUMENT STERILIZING TRAY

Material : Polypropylene

This autoclavable trays have cover which fits into groove of trays to minimize spillage of material. Smaller size can be used to sterilize small instruments whereas the larger size allows pipettes to be fully immersed in disinfectant prior to sterilizing.

| Product Code | Particulars | HSN Code |
|--------------|--------------|----------|
| 81736 | 220x150x70mm | 39269099 |
| 81737 | 450x150x70mm | 39269099 |



ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

DRAINING RACK

Material : Polypropylene

Poly;ab Draining Rack provides an excellent option for post wash draining of tubes & bottles in the laboratory. The back plate has holes into which pegs are tightened with the help of screw-nuts. Bottom of the rack is provided with a tray which collect all the water drained down which can further be collected in a container via drain tube. This rack does not, stain or rust.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 81741 | 20 Pegs | 39269099 |
| 81742 | 56 Pegs | 39269099 |



SINCE 1966

MAGENTA BOX

Material : Polycarbonate/Polyprpylene

Magenta Box is very often used in tissue culture & agricultural research experiments. main body of this box is made of tough & transparent Polycarbonate which provides excellent view of the contents inside while the cover is made of Polypropylene. This autoclavable box is square in shape giving itself a space saving feature.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 81801 | Magenta Box | 39231090 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS



STORAGE BOXES

Material : Polypropylene

These Storage Boxes, made of polypropylene, prove to be an excellent option for storage of different powdered reagents & other material in the laboratory. Square shape of these boxes help in saving shelf space. Available in three different sizes, these boxes have polypropylene caps, are autoclavable, provide good chemical resistance & heavy duty services.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 81900 | 250ml | 39231090 |
| 81901 | 500ml | 39231090 |
| 81902 | 1000ml | 39231090 |



MULTI PURPOSE RECTANGULAR STORAGE BOX

Material : PP Material Conforming to US FDA 21CFR.

These multipurpose storage boxes are manufactured from polypropylene for storing and transporting all types of products including food stuff.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 81911 | Storage Box | 39321090 |

SINCE 1966



SIMPLE CELL POT

Material : Polypropylene

Simple Cell Pot, as the name suggests, is used for making Simple Cells & Voltmeters etc. Moulded in Polypropylene, this Pot has very good contact clarity and excellent chemical resistance. Being in Polypropylene, this Pot is much durable that the one used conventionally.

| Product Code | Particulars | HSN Code |
|--------------|-----------------|----------|
| 82101 | Simple Cell Pot | 39269099 |



ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

LECLANCHE CELL POT

Material : Polypropylene

Blown in Polypropylene, this Pot is used for making Leclanche Cell. This Pot too has very good contact clarity & excellent chemical resistance.

| Product Code | Particulars | HSN Code |
|--------------|--------------------|----------|
| 82102 | Leclanche Cell Pot | 39269099 |



SAMPLE CONTAINER (PRESS & FIT TYPE)

Material : PP Material Conforming to US FDA 21CFR.

Made of polypropylene, these disposable containers are used for collecting samples of sputum or urine etc. Graduated on the outer surface, these containers have cap that work in a press & fit manner & make these containers leakproof. These disposable containers are individually packed and are E.T.O. sterilized.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 82103 | 50 ml | 39233090 |
| 82104 | 25ml | 39233090 |



SINCE 1966

MEDICINE CUP

Material : PP Material Conforming to US FDA 21CFR.

These transparent polypropylene cups can be used to measure different volumes of medicines to be administered. The inbuilt graduations make measuring very easy.

| Product Code | Particulars | HSN Code |
|--------------|-------------|----------|
| 83101 | 25ml | 39269099 |



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

PHYSICAL PROPERTIES OF POLYLAB LABWARE

| RESIN | Max use Temp ^o C | Bittleness Temp ^o C | Transparency | Sterilization | | | | | | Flex-ibility | Water Absorp-tion% |
|-------------|-----------------------------|--------------------------------|--------------|---------------|-----|----------|------------|---------------|------------------|--------------|--------------------|
| | | | | Auto-claving | Gas | Dry-Heat | Radia-tion | Disin-fectant | Specific Gravity | | |
| HDPE | 120 | -100 | Translucent | No | Yes | No | Yes | Yes | 0.95 | Rigid | 20.01 |
| LDPE | 80 | -100 | Translucent | No | Yes | No | Yes | Yes | 0.92 | Excel | <0.01 |
| PC | 135 | -135 | Clear | Yes+ | Yes | No | No | Yes | 1.20 | Rigid | 0.35 |
| PP | 135 | 0 | Translucent | Yes | Yes | No | No | Yes | 0.90 | Rigid | <0.02 |
| PS | 90 | 20 | Clear | No | Yes | No | Yes | Some | 1.05 | Rigid | 0.05 |

Sterilization

- Autoclaving (121^oC, 15 psig for 20 minutes) Clean and rinse item with distilled water before autoclaving. Certain chemicals which have no appreciable effect on resins at room temperature may cause deterioration at autoclaving temperatures.
 - ALWAYS COMPLETELY DISENGAGE THREADS BEFORE AUTOCLAVING.
 - Gas - Ethylene oxide formaldehyde.
 - Dry Heat - 160^oC for 120 minutes
 - Disinfectant - Benzalkonium chloride, formalin, ethanol, etc.
 - Radition - gamma irradiated at 2.5 Mard with unstabilized plastic.
- Sterilizing reduces mechanical strength. Do not use PC vessels for vaccum application if they have been Autoclaved.



ASIAN SCIENTIFIC INSTRUMENTS
SUPER MARKET FOR LAB ITEMS

CHEMICAL RESISTANCE FOR POLYLAB LABWARE

| CLASSES OF SUBSTANCES AT ROOM TEMPERATURE | HDPE | LDPE | PC | PP | PS |
|---|--------|--------|--------|--------|--------|
| ACIDS, DILUTE OR WEAK | Green | Green | Green | Green | Green |
| ACIDS, STRONG AND CONCENTRATED | Green | Green | Red | Green | Orange |
| ALCOHOLS, ALIPHATIC | Green | Green | Green | Green | Green |
| ALDEHYDES | Blue | Blue | Orange | Blue | Red |
| BASES | Green | Green | Red | Green | Green |
| ESTERS | Blue | Blue | Red | Blue | Red |
| HYDROCARBONS, ALIPHATIC | Blue | Orange | Orange | Blue | Red |
| HYDROCARBONS, AROMATIC | Blue | Orange | Red | Orange | Red |
| HYDROCARBONS, HALOGENATED | Orange | Red | Red | Orange | Red |
| KETONES | Blue | Blue | Red | Blue | Red |
| OXIDIZING AGENTS, STRONG | Orange | Orange | Red | Orange | Red |

- 30 Days of constant exposure cause no damage. Plastics may even tolerate for years.
- Little or no damage after 30 days of constant exposure to the reagents.
- Some effect after 7 days of constant exposure to the reagent like crazing, cracking, loss of strength or discoloration.
- Not recommended for continuous use, immediate damage may occur.



✉ info@asianhyd.com

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

ABBREVIATIONS

| | |
|-------------|---|
| PP | Polypropylene |
| PC | Polycarbonate |
| HDPE | High Density Polyethylene |
| LDPE | Low Density Polyethylene |
| PS | Polystyrene |
| ABS | Acrylo Nitrile Butadiene Styrene |
| PTFE | Poly Tetra Flouro Ethylene |
| PMMA | Poly Methyl Methacrylate |
| TPX | Poly Methyl Pentene |

TERMS & CONDITIONS

- **G.S.T** will be charged extra as applicable at the time of billing.
- **The Prices** are subject to alteration without any prior notice.
- **The Price** for items of special requirements will be quoted on request.
- **Packing and forwarding** are extra as applicable.
- **Transportation Charges** are on to Pay Basis.
- **Our responsibility** ceases once we deliver the goods to the carriers.
- **Goods** will be dispatched in standard packing or multiples thereof.
- **All disputes** shall be within the jurisdiction of the competent courts in Delhi
- **Capacity measures** are for laboratory use only and not for commercial use.
- **This cancels** all our previous Price List.
- **The company** reserves the right to change colour & modify the design of the product.
- **This Price List** is effective from 1st October 2023

✉ info@asianhyd.com





ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

ASIAN SCIENTIFIC INSTRUMENTS

Asian Estate : 4-1-10, Tilak Road Opposite Posnett Bhavan Hyderabad, Telangana-500 001
Ph.: 93930-01939, 90522-23905, 98480-31317, 93472-45096

✉ info@asianhyd.com