

DIFFERENTIAL BLOOD CELL COUNTER

“ISO 9001:2008 & ISO 13485:2003”
CERTIFIED COMPANY
“CE MARKED PRODUCTS”



Model No : DC - 22

Mains Input	: 230V/50Hz 110V/60Hz A.C. Supply
Body	: ABS Body (Light Weight)
LCD Display	: 16 Character With Back Light
Weight	: 0.6 Kg.
Overall Size	: 155.0 (L) x 205 (W) x 91 (H)mm
Cell Counting Key	: 8 KEY
Total Key	: 12 KEY
Keypad	: Tactile switch with raised keytops
Timer & Stop Upto	999 Minutes