




DIGITAL BALANCES



KRT Series are entry-level precision, compact Balance with rich features and suitable for Laboratory & academic applications. The performance that users have come to expect from **Kerov**

MODEL	KRT 200	KRT 300	KRT 600	KRT 1000H	KRT 1000S	KRT 300S	KRT 1000	KRT 2000	KRT 3000
Capacity (g)	210	310	610	1000	1100	310	1000	2100	3100
Readability (g)	0.01			0.02	0.005	0.1			
Repeatability (g) ±	0.01			0.02	0.005	0.1			
Linearity (g) ±	2 count								
Tare Range	100%								
Stabilization Time	3 second								
Operating Temp	10-40 C								
Pan Size (mm) Approx.	Ø116						Ø150		
Power Requirements	AC 220V 50Hz								
Windshield	Standard						Optional		

**You can contact us for any Capacity and Readability.*