


Print



Model	KRT200	KRT300	KRT600	KRT1000	KRT2000	KRT3000	KRT1000H
Capacity (g)	200	300	600	1000	2000	3000	1000
Readability (g)	0.01						0.5
Repeatability (g)	+/- 0.01						+/- 0.05
Linearity (g)	+/- 0.02						+/- 0.05
Tare Range	within range						
Stabilization Time	3 seconds						
Operating Temp.	10 - 45 C						
Pan Size (mm)	116				150	116	
Dimensions (mm)	260 x 185 x 75						
Power Requirement	AC 220V 50 Hz.						