

SISCOTM

PREMIUM +

An ISO 9001-2008 & CE Certified Company



MUFFLE FURNACE



SISCOTM

INSTRUMENTS

8UA, Jawahar Nagar, DELHI-110007
Email id: siscoinstruments@gmail.com



MUFFLE FURNACE

Specifications :

- The outer wall made of Mild steel sheet powder coated.
- High grade ceramic wool insulation to minimize thermal loss.
- Design Temp. range up to to 1200°C Working temp 1100°C controller by a Microprocessor PID based Digital Temperature Indicator cum Controller.
- Heating Element made of Kanthal wire special high temperature alloy to ensure uniform distribution of heat.
- To work on 220/440 volt, 50Hz. Single or Three Phase Supply.

	Inner chamber size	Model	Working Temperature upto	Rating (K.W.) Approx
A	225x100x100 mm (9"x4"x4")	Pharmacy	900 °C	2.5
B	250x125x125 mm (10"x5"x5")	Industrial	1200 °C	2.5
C	300x150x150 mm (12"x6"x6")	Industrial	1200 °C	4.0

***Other Sizes also available on request.

Optional : a) Profiler controller with Ramp and sock facility 16 Steps



Note : SISCO Premium GMP Model with complete inner and outer Stainless Steel also available.

Authorized Distributor:-