


HEATING mantle



**“ WHERE
PRECISION MEETS
PERFORMANCE ”**

In scientific research and chemical laboratories, accurate temperature control is essential. Introducing the SISCO Instruments Heating Mantle—your reliable partner in achieving consistent, safe, and efficient heating solutions for laboratory glassware. Designed with precision, built for durability.



**Best in
Accurate
temperature**

**START YOUR RESEARCH
WITH THE PERFECT ONE**

HEATING mantle

Why choose **SISCO** heating Mantle?



Heating Mantles offer safe, uniform, and efficient heating for all your lab needs. With durable construction, glass wool insulation, and precise temperature control, they ensure reliable performance every time. Trusted by professionals for quality and consistency.

FEATURES

- Uniform heating with energy efficiency
- Built-in energy regulator for temperature control
- Strong insulation with glass wool for safety
- Robust design with glass sleeving reinforcement
- Available in various capacities (250 ml to 10,000 ml)
- Easy to operate and maintain
- Ideal for long-duration lab use



Sisco Instruments Heating Mantles are designed for safe, uniform, and efficient heating of laboratory glassware. With glass sleeving and quality insulation, they provide consistent performance and protection during use.



Safe



Efficient



Uniform heat



wide Range

What we offer

- High-Performance Heating Mantles:** Built with durable glass sleeving and insulated with premium glass wool, our mantles ensure safe, uniform, and efficient heating for all lab applications.
- Precision Temperature Control:** Each unit comes with a built-in energy regulator, allowing smooth and accurate control over heating, reducing the risk of overheating or inconsistencies.
- Versatile Range of Models:** Available in various capacities from 250 ml to 10,000 ml and power ratings up to 1000W, we offer solutions to meet every lab's unique needs.



APPLICATIONS

- Chemistry and pharmaceutical labs
- Industrial R&D departments
- Biological and chemical reaction experiments
- Evaporation and distillation processes
- Educational institutions for practical experiments




HEATING INDICATOR

ENERGY REGULATOR



HEATING mantle

Choose right one



Choosing the right heating mantle depends on your lab's specific heating requirements. Whether you're working with small samples or large volumes, having the correct capacity and wattage ensures efficient and uniform heating.

Our range of models offers flexibility for various applications, with each unit built for accuracy, durability, and safety. Explore the technical specifications below to find the model that best fits your needs.

Available models

POWER	CAPACITY
150 W	250 ml
200 w	500 ml
300 w	1000 ml
450 w	2000 ml
500 w	3000 ml
700 w	5000 ml
1000 w	10000 ml
1800 w	20000 ml



Disclaimer :

The images in this brochure are for reference purposes only and may differ from the actual product or service.

For inquiries, please contact us at:

siscoinstruments@gmail.com