

# SISCO<sup>TM</sup> PREMIUM +

An ISO 9001-2008 & CE Certified Company



## B.O.D. INCUBATOR



SISCO<sup>TM</sup>  
INSTRUMENTS

8UA, Jawahar Nagar, DELHI-110007  
Email id: [siscoinstruments@gmail.com](mailto:siscoinstruments@gmail.com)



# B.O.D. INCUBATOR

## SPECIFICATIONS:

- Double Walled Inside made of Stainless Steel
- Outside made of M.S. Sheet painted in epoxy powder coating.
- The gap between outer and inner Wall is filled high grade thermal insulation.
- For Temperature uniformity Air circulation 2 Nos. fitted inside.
- Digital Voltmeter, castor wheels for moment.
- Temp. range +5°C to 50°C By Digital PID Temperature Indicator cum Controller with an Temperature Accuracy  $\pm 1^\circ\text{C}$ .
- To work on 220/230 Volt AC. Supply.

	Inner Chamber Size in mm (HxWxD)	Capacity cu.ft.	Capacity Liters.	Connected Kw
1	455x410x610	4.1	112	1
2	505x415x830	6.1	171	1.5
3	565x565x865	10.1	280	2

\*\*\*Other Sizes also available on request.

### Optional :

- a) Co2 Air mixer Nozzle is Provide.
- b) Internal illumination with 3 fluorescent tube.
- c) Automatic Voltage Stab with 0-24 timer.
- d) PLC controller with 21 CFR software and data recorder

### B.O.D. Incubator (Microprocessor Based)



### Note :

SISCO Premium GMP Model with complete inner and outer Stainless Steel also available.

**Authorized Distributor:-**



**ASIAN SCIENTIFIC INSTRUMENTS**

**SUPER MARKET FOR LAB ITEMS**

Asian Estate, 4-1-10, Opp. Posnett Bhavan, Tilak Road, Hyderabad-500 001.

☎ : 90522 23905, 93930 01939, 92810 93814, 93999 62773, 95502 60447

☎ : 040-2475 1902 / 2475 4607

✉ : info@asianhyd.com 🌐 www.asiansci.com