

SISCOTM PREMIUM +

An ISO 9001-2008 & CE Certified Company



HOT PLATE



SISCOTM
INSTRUMENTS

8UA, Jawahar Nagar, DELHI-110007

Email id: siscoinstruments@gmail.com



HOT PLATE

SPECIFICATIONS :-

- Outer body made of mild steel powder coated, Top made of Cast iron.
- Temperature range up to 350°C With Kanthal Element are placed underneath for uniform heating.
- To work on 220/230 Volt 50Hz AC Supply.
- Supplied with Energy Regulator (Temp range 50-300°C)
- MS Powder coated working Top Platform

	SIZE	LOAD (kw)
A	20 cm Round/ Single Stage (8")	1
B	300x250 mm (12"x10")	1
C	400x250mm (16"x10")	1.5
D	450x300mm (18"x12")	2
E	600x450mm (24"x18")	3

***Other Sizes also available on request.



Optional : a) Microprocessor based PID Digital Temp. Indicator cum Controller (Instead of Energy Regulator)

Note : SISCO Premium GMP Model with complete inner and outer Stainless Steel also available.



HOT PLATE

SPECIFICATIONS :-

- Outer body made of mild steel powder coated, Top made of Stainless steel SS304.
- Temperature ranges up to 350°C With Kanthal Element are placed underneath for uniform heating.
- To work on 220/230 Volt 50Hz AC Supply.
- Supplied with Energy Regulator (Temp range 50-300°C)
- Stainless steel working Top Platform

	SIZE
A	20 cm Round/ Single Stage (8")
B	20 cm Round/ Double Stage



Note :

SISCO Premium GMP Model with complete inner and outer Stainless Steel also available.

Authorized Distributor:-



SINCE 1966

ASIAN SCIENTIFIC INSTRUMENTS

SUPER MARKET FOR LAB ITEMS

Asian Estate, 4-1-10, Opp. Posnett Bhavan, Tilak Road, Hyderabad-500 001.

☎ : 90522 23905, 93930 01939, 92810 93814, 93999 62773, 95502 60447

☎ : 040-2475 1902 / 2475 4607

✉ : info@asianhyd.com 🌐 www.asiansci.com